REF. 2138 289 97 I

FR	NOTICE D'UTILISATION – IMPLANTS SCLERAUX POUR LE DECOLLEMENT DE RETINE	2
EN	INSTRUCTION FOR USE – SCLERAL BUCKLING IMPLANTS FOR RETINAL DETACHMENT	3
DE	GEBRAUCHSANWEISUNG – IMPLANTAT ZUR SKLERALEN BUCKELUNG BEI NETZHAUTABLÖSUNG	4
IT	ISTRUZIONI PER L'USO – IMPIANTISCLERALI PER IL DISTACCO DELLA RETINA	5
ES	INSTRUCCIONES DE USO – IMPLANTES ESCLERÓTICOS PARA EL DESPRENDIMIENTO DE RETINA	6
PT	INSTRUÇÕES DE USO – IMPLANTES ESCLERÓTICOS PARA O DESCOLAMENTO DA RETINA	7
NL	GEBRUIKSAANWIJZING – SCLERALE IMPLANTATEN VOOR HET LOSWEKEN VAN HET NETVLIES	8
DA	BRUGSVEJLEDNING – IMPLANTATER TIL SCLERAL BUCKLING VED RETINALØSNING	9
FI	KÄYTTÖOHJEET – IMPLANTAT FÖR SKLERAL UTBUKTNING VID NÄT-HINNEAVLOSSNING	10
SE	BRUKSANVISNING – IMPLANTAT FÖR SKLERAL UTBUKTNING VID NÄT-HINNEAVLOSSNING	11
TR	KULLANIM – RETİNAL DEKOLMAN AMAÇLI SKLERAL ÇÖKERTME İMPLANTLARI	12
ZH	使用说明 – 视网膜脱离治疗用巩膜扣带植入物	13
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - ΕΜΦΥΤΕΥΜΑΤΑ ΣΚΛΗΡΙΑΙΑΣ ΕΠΑΝΑΣΥΓΚΟΛΛΗΣΗΣ ΓΙΑ ΤΗΝ ΑΠΟΚΟΛΜΗΣΗ ΤΟΥ ΑΜΦΙΒΛΗΣΤΡΟΕΙΔΟΥΣ	14
PL	INSTRUKCJA UŻYTKOWANIA – IMPLANTY DO KLAMROWANIA TWARDÓWKI W PRZYPADKU ODWARSTWIENIA SIATKÓWKI	15
SL	NAVODILA ZA UPORABO – VSADKI ZA UPOGIB BELOČNICE PRI ODPSTOPU MREŽNICE	16
SK	INSTRUCTION FOR USE – SCLERAL BUCKLING IMPLANTS FOR RETINAL DETACHMENT	17

Figures 18

Symbols 19

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EN

Instructions for use
Date of revision of the instructions for use : 06/2016

SCLERAL BUCKLING IMPLANTS FOR RETINAL DETACHMENT

Description:

The SCLERAL BUCKLING IMPLANTS are available in different shapes, each further divided into different sizes:

- Bands
- Fixation Sleeves : designed to help the fixation of a band during a complete encircling procedure. Available in different sizes depending on the band model.
- Strips : designed for placement under bands. Available in different models depending on the width and height of the indentation required.
- Tires : designed for placement under bands for greater indentation, in cases of multiple and wide tearing or of vitreoretinal proliferation.
- Cords
- Sponges, with the following types of sections:
 - round
 - oval
 - grooved
 - half-round

These implants are available in the following different materials:

- Solid Silicone
- Silicone sponge
- Expanded polytetrafluoroethylene (ePTFE)

Packaging:

SCLERAL BUCKLING IMPLANTS are supplied sterile. They are sterilized in their final double packaging to make them easier to handle in aseptic conditions. The silicone implants are sterilized by gamma irradiation (minimal dose : 25 kGy). A red indicator shows that the product followed a sterilization cycle validated by FCI. Do not use the product if the indicator is not red.

The ePTFE implants are sterilized using ethylene oxide. A green indicator shows that the product followed a sterilization cycle validated by FCI. Do not use the product if the indicator is not green.

The method of sterilization of each implant is indicated on the packaging.

Indications:

The SCLERAL BUCKLING IMPLANTS are used for indentation of the sclera during retinal detachment surgery. They enable a re-application of the choroid to the retina. They are produced from non resorbable materials and provide permanent indentation. The different materials available (Solid silicone, silicone sponge, ePTFE) offer different mechanical properties of elongation, compressibility and elasticity, thus allowing the surgeon to obtain the right degree of indentation.

Contraindication:

Pre-existing physiological characteristics or pathologies that can be exacerbated by the implantation of SCLERAL BUCKLING IMPLANTS are the following :

- major scleral thinning,
- conjunctivitis,
- septic risk.

Complications:

As in all types of surgery, there are risks linked to the material used or to developments in the initial pathology. Potential complications associated with retinal surgery detachment or scleral indentation can be the following :

Intraoperative complications :

- appearance of a fishmouth

Post-operative :

- detachment of the choroid
- exteriorization of the implant
- scleral erosion by the implant
- edema of the corneal epithelium
- angle-closure glaucoma
- secondary retinal detachment
- ischemia of the anterior segment
- infection, endophthalmitis
- strabismus
- cystoid macular edema

Unexpected side effects and complications related to the SCLERAL BUCKLING IMPLANTS must be reported to FCI.

Warnings:

SCLERAL BUCKLING IMPLANTS must be removed from their packaging and handled in aseptic conditions. Before use, the individual pack preserving the product's sterility should be checked to make sure it is intact. SCLERAL BUCKLING IMPLANTS are single-use products and must not be re-sterilized. Reuse and/or re-sterilization may compromise device performance, which could cause serious harm to the patient's health and safety. They should be stored at room temperature and must not be used after the expiration date shown on the package.

The use of Watzke forceps S5.3100 is recommended for the installation of fixation sleeves.

The guaranteed maximum residual rate of ethylene oxide at the time the device was placed on the market was 5.3 µg per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

An implant card is available in the product's packaging. It is intended to be provided to the patient and allows traceability of the implanted product.

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