REF. 2138 299 34 L

FR	NOTICE D'UTILISATION - LENTILLES À USAGE UNIQUE POUR VITRECTOMIE	2
EN	INSTRUCTION FOR USE - DISPOSABLE VITRECTOMY LENSES	3
DE	GEBRAUCHSANWEISUNG - EINWEGLINSEN FÜR DIE VITREKTOMIE	4
IT	ISTRUZIONI PER L'USO - LENTI OCULARI USA E GETTA PER VITRECTOMIA	5
ES	INSTRUCCIONES DE USO - LENTES DE USO ÚNICO PARA VITRECTOMÍA	6
PT	INSTRUÇÕES DE USO - LENTES DE USO ÚNICO PARA VITRECTOMIA	7
NL	GEBRUIKSAANWIJZING - WEGWERPLENZEN VOOR VITRECTOMIE	8
DA	BRUGSVEJLEDNING - VITREKTOMILINSE TIL ENGANGSBRUG	9
FI	KÄYTTÖOHJEET - KERTAKÄYTTÖINEN VITREKTOMIALINSSI	10
SE	BRUKSANVISNING - VITREKTOMILINS FÖR ENGÅNGSBRUK	11
TR	KULLANIM - TEK KULLANIMLIK VITREKTOMİ LENSLERİ	12
ZH	使用说明 - 玻璃体切割术用一次性透镜	13
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - ΦΑΚΟΙ ΥΑΛΟΕΙΔΕΚΤΟΜΗΣ ΜΙΑΣ ΧΡΗΣΗΣ	14
PL	INSTRUKCJA UŻYTKOWANIA – JEDNORAZOWE SOCZEWKI DO WITREKTOMII	15
SL	NAVODILA ZA UPORABO - LEČE ZA VITREKTOMIJO ZA ENKRATNO UPORABO	16
SK	NÁVOD NA POUŽITIE - JEDNORAZOVÉ VITREKTOMICKÉ ŠOŠOVKY	17

Symbols.....18

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DISPOSABLE VITRECTOMY LENSES

Description :

The DISPOSABLE VITRECTOMY LENSES are made of optical grade silicone.

Packaging :

The DISPOSABLE VITRECTOMY LENSES are supplied sterile (sterilized by ethylene oxide). A green indicator shows that the product followed a sterilization cycle validated by FCI. Do not use the product if the indicator is not green.

For ease of handling in aseptic conditions, the DISPOSABLE VITRECTOMY LENSES are presented in a double packaging.


Indications :

The DISPOSABLE VITRECTOMY LENSES are used in the context of vitrectomy surgeries.

The lenses, used with a microscope, facilitate visualisation of the retina by cancelling the corneal dioptré.

An edge, which has a frosted ring on its upper side in order to avoid edge glare effect, allows users to position and move the lens on the cornea with forceps.

The top of this lip is left unpolished in order to distinguish it from the polished side that is placed on the cornea.

Warning and contra-indications:

Any secondary effects or unexpected complications possibly linked to the use of the DISPOSABLE VITRECTOMY LENSES must be reported to FCI.

Recommendations for use :

To be used with a microscope.

The surfaces are polished and care must be taken to ensure that they are not damaged or dirtied.

Before applying the lens against the cornea, lubricate its concave side with a drop of ophthalmic gel.

Precautions for use :

The DISPOSABLE VITRECTOMY LENS should be removed from its packaging and handled in aseptic conditions. Check the individual sterilisation protection prior to use. The DISPOSABLE VITRECTOMY LENSES are single use products and must not be re-sterilised. Reuse and/or reesterilization may compromise device performance, which could cause serious harm to the patient's health and safety.

It should be kept at room temperature and must not be used after the expiration date stated on the packaging.

The guaranteed maximum residual rate of ethylene oxide at the time the device was placed on the market was 5.3 µg per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

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