REF. 2138 299 77 H

FR	NOTICE D'UTILISATION - BILLES BIOCERAMIQUES OVOIDES	2
EN	INSTRUCTION FOR USE - BIOCERAMIC EGG-SHAPED ORBITAL IMPLANTS	3
DE	GEBRAUCHSANWEISUNG - EIFÖRMIGES ORBITAIMPLANTAT AUS BIOKERAMIK	4
IT	ISTRUZIONI PER L'USO - SFERE IN BIOCERAMICA OVOIDALI	5
ES	INSTRUCCIONES DE USO - BOLAS BIOCERÁMICAS OVOIDES	6
PT	INSTRUÇÕES DE USO - ESFERAS BIOCERÁMICAS OVÓIDES	7
NL	GEBRUIKSAANWIJZING - BIOKERAMISCHE EIVORMIGE IMPLANTATEN	8
DA	BRUGSVEJLEDNING - ÆGGEFORMEDE BIOKERAMISKE ØJEIMPLANTATER	9
FI	KÄYTTÖOHJEET - BIOKERAAMISET MUNANMUOTOISET ORBITAIMPLANTIT	10
SE	BRUKSANVISNING - ÄGGFORMAT ORBITAIMPLANTAT AV BIOKERAMIK	11
TR	KULLANIM - YUMURTA ŞEKLİNDE BİYOSERAMİK ORBİTAL İMPLANTLAR	12
ZH	使用说明 - 义眼台 (椭圆型生物陶瓷)	13
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - ΩΟΕΙΔΗ ΕΝΔΟΚΟΓΧΙΚΑ ΕΜΦΥΤΕΥΜΑΤΑ ΑΠΟ ΒΙΟΚΕΡΑΜΙΚΟ ΥΛΙΚΟ	15
PL	INSTRUKCJA UŻYTKOWANIA – BIOCERAMICZNE OWALNE IMPLANTY OCZODOŁOWE	16
SL	NAVODILA ZA UPORABO – BIOKERAMIČNI ORBITALNI VSADKI V OBLIKI JAJCA	17
SK	ÁVOD NA POUŽITIE - BIOKERAMICKÉ ORBITÁLNE IMPLANTÁTY V TVARE VAJÍČKA	18

Figures	19
Symbols	20

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BIOCERAMIC EGG-SHAPED ORBITAL IMPLANTS

Description and packaging:

BIOCERAMIC EGG-SHAPED ORBITAL IMPLANTS are intraorbital alumina implants with 99 % guaranteed purity. Their average porosity is 250-550 µm. They are available in several diameters : 18-20-21 mm (S6.5330, S6.5340, S6.5350).

The BIOCERAMIC EGG-SHAPED ORBITAL IMPLANTS are supplied sterile. They are sterilized using ethylene oxide. A green indicator shows that the product followed a sterilization cycle validated by FCI.

BIOCERAMIC ORBITAL IMPLANTS are sterilized in their final double packaging, to make them easier to handle in aseptic conditions.

Indications :

The BIOCERAMIC EGG-SHAPED ORBITAL IMPLANTS are designed to fill in the orbital cavity following enucleation, evisceration or during secondary implantation.

Mode of action :

The porosity of the BIOCERAMIC EGG-SHAPED ORBITAL IMPLANT facilitates its colonization by fibrovascular tissue, which offers the following advantages:

- reduced risk of infection,
- reduced risk of implant extrusion,

Contraindications :

Care should be given to irradiated orbital cavities where the capacity for vascularization and colonization can be reduced (Especially in cases of irradiation of the orbital cavity and in diabetic patients).

Adverse effects :

As in any type of surgery, there are risks linked to the material or to developments of the initial pathology. Potential complications associated with the implantation of a BIOCERAMIC EGG-SHAPED ORBITAL IMPLANT include, but are not limited to the following :

- Conjunctival dehiscence
- Local infection or irritation
- Exposure
- Discomfort
- Ptosis
- Laxity of eyelids
- Enophthalmos
- Hypo-ophthalmia

Unexpected side effects and complications related to the BIOCERAMIC EGG-SHAPED ORBITAL IMPLANTS must be reported to FCI.

Advice for use :

The choice of the diameter is an important parameter for the final result. Given the retractile effect of orbital tissues, and to avoid dehiscence around the scar, the BIOCERAMIC EGG-SHAPED ORBITAL IMPLANT diameter has to be adapted to the patient's anatomy.

The use of S6.3060 sizers is advisable. The proper diameter should allow the wound to be closed without tension. The choice of the diameter must take into account for the thickness of any material used to wrap the implant.

The use of sphere introducer S6.3050 is recommended for the install of the BIOCERAMIC EGG-SHAPED ORBITAL IMPLANT.

The enucleation or the evisceration technique is conventional

BIOCERAMIC EGG-SHAPED ORBITAL IMPLANTS can be implanted uncovered, particularly in cases of evisceration.

In cases of evisceration, it is usually not advisable to retain the cornea, particularly if there is history of corneal injury.

In cases of enucleation : it is possible to cover the BIOCERAMIC EGG-SHAPED ORBITAL IMPLANT with either absorbable mesh or autogenous material (sclera, fascialata, temporal aponeurosis).

Precautions for use :

Remove the BIOCERAMIC EGG-SHAPED ORBITAL IMPLANT from its packaging in aseptic conditions. Before use, verify that the individual pack preserving the product's sterility is intact. BIOCERAMIC EGG-SHAPED ORBITAL IMPLANT are single-use products that should not be resterilized. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety. Do not use the product if the indicator is not green, as it could endanger the patient. They should be stored at room temperature and should not be used after the expiration date shown on the package.

The drilling of the sphere is not recommended as the process may damage it.

The guaranteed maximum residual rate of ethylene oxide at the time the device was placed on the market was 5.3 µg per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

An implant card is available in the product's packaging. It is intended to be provided to the patient and allows traceability of the implanted product.

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