REF. 2138 302 90 D

FR	NOTICE D'UTILISATION - VERRE DE CONTACT TROIS MIROIRS A USAGE UNIQUE	2
EN	INSTRUCTION FOR USE - SINGLE-USE THREE-MIRROR CONTACT LENS	3
DE	GEBRAUCHSANWEISUNG - DREISPIEGELKONTAKTGLAS ZUM EINMALGEBRAUCH	4
IT	ISTRUZIONI PER L'USO - LENTE A CONTATTO TRE SPECCHI USA E GETTA	5
ES	INSTRUCCIONES DE USO - VIDRIO DE CONTACTO TRES ESPEJOS PARA USO ÚNICO	6
PT	INSTRUÇÕES DE USO – VIDRO DE CONTACTO TRÊS ESPELHOS PARA USO ÚNICO	7
NL	GEBRUIKSAANWIJZING - DRIE SPIEGELS CONTACTLENS VOOR EENMALIG GEBRUIK	8
DA	BRUGSVEJLEDNING - TRESPEJLSKONTAKTLINSE TIL ENGANGSBRUG	9
FI	KÄYTTÖOHJEET - KERTAKÄYTTÖINEN KOLMEN PEILIN LINSSI	10
SE	BRUKSANVISNING - TRESPEGELKONTAKTGLAS FÖR ENGÅNGSBRUK	11
TR	KULLANIM - TEK KULLANIMLIK ÜÇ AYNALI KONTAKT LENS	12
ZH	使用说明 - 一次性三面镜	13
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - ΤΡΙΚΑΤΟΠΤΡΙΚΟΣ ΦΑΚΟΣ ΕΠΑΦΗΣ ΜΙΑΣ ΧΡΗΣΗΣ	14
PL	INSTRUKCJA UŻYTKOWANIA – JEDNORAZOWA SOCZEWKA KONTAKTOWA Z TRZEMA LUSTERKAMI	15
SL	NAVODILA ZA UPORABO – TROZRCALNA KONTAKTNA LEČA ZA ENKRATNO UPORABO	16
SK	NÁVOD NA POUŽITIE - JEDNORAZOVÉ TROJZRKADLOVÉ KONTAKTNÉ ŠOŠOVKY	17

Symbols..... 18

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SINGLE-USE THREE-MIRROR CONTACT LENS

Description :


THE SINGLE-USE THREE -MIRROR CONTACT LENS comprises four viewing zones.

Zone 1: Lens for observing the posterior pole up to 30°.

Zone 2: 75° inclined mirror for observing the fundus above the 30° mark.

Zone 3: 67° inclined mirror for observing the peripheral retina.

Zone 4: 59° inclined mirror for observing the vitreous body and the areas of the fundus alongside the ora serrata, as well as for gonioscopic examinations.

Packaging :

THE SINGLE-USE THREE -MIRROR CONTACT LENS is sterilized using ethylene oxide and packaged in double packaging. A green indicator shows that the product followed a sterilization cycle validated by FCI. Do not use the product if the indicator is not green.

Indications :

THE SINGLE-USE THREE -MIRROR CONTACT LENS is used for examining the fundus of the eye and the iridocorneal angle.

Complications :

Unexpected side effects and complications that can be attributed to the SINGLE-USE THREE -MIRROR CONTACT LENS must be reported to FCI. There is no contraindication to use the contact lens. In case of ocular disease or corneal injuries, the practitioner must take additional precautions according to the clinical case.

THE SINGLE-USE THREE -MIRROR CONTACT LENS is made out of acrylic must never be used for laser treatment.

Intructions for use :

When used in conjunction with a slit lamp or a microscope during an operation, this contact lens enables the retina to be viewed by eliminating the corneal diopter. You are advised to ensure that the pupils are considerably dilated at the time of use.

The surfaces of the SINGLE-USE CONTACT LENS are perfectly smooth and care must be taken not to damage or dirty them when handling.

Remove the lens from its packaging under aseptic conditions. Before applying against the cornea, lubricate the concave side of the lens with a drop of ophthalmic gel.

Warnings :

The CONTACT LENS must be stored at room temperature and must not be used after the expiration date indicated on the packaging. Check that the individual sterile protector is intact before use. THE SINGLE-USE THREE -MIRROR CONTACT LENS IS FOR SINGLE USE. Reuse and/or reesterilization may compromise device performance, which could cause serious harm to the patient's health and safety.

The guaranteed maximum residual rate of ethylene oxide at the time the device was placed on the market was 5.3 µg per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

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