REF. 2138 307 83 H

FR	NOTICE D'UTILISATION - PAINLESS PLUG	2
EN	INSTRUCTION FOR USE - PAINLESS PLUG / SNUG PLUG	3
DE	GEBRAUCHSANWEISUNG - PAINLESS PLUG	4
IT	ISTRUZIONI PER L'USO - PAINLESS PLUG	5
ES	INSTRUCCIONES DE USO - PAINLESS PLUG	6
PT	INSTRUÇÕES DE USO - PAINLESS PLUG	7
NL	GEBRUIKSAANWIJZING - PAINLESS PLUG	8
DA	BRUGSVEJLEDNING - SMERTEFRI PROP	9
FI	KÄYTTÖOHJEET - KIVUTON TULPPA	10
SE	BRUKSANVISNING - SMÄRTFRI PLUGG	11
TR	KULLANIM - AĞRISIZ TIKAÇ	12
ZH	使用说明 - 泪点塞(painless 式)	13
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - PAINLESS PLUG / SNUG PLUG	15
PL	INSTRUKCJA UŻYTKOWANIA – ZATYCZKA PAINLESS PLUG / SNUG PLUG	16
SL	NAVODILA ZA UPORABO – BREZBOLEČINSKI ČEP/TESNI ČEP	17
SK	NÁVOD NA POUŽITIE - BEZBOLESTNÁ ZÁTKA/PRILIEHAVÁ ZÁTKA	18

Figures 19

Symbols 20

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PAINLESS PLUG / SNUG PLUG

Description and packaging:

The PAINLESS PLUG / SNUG PLUG is a punctal plug made of silicone and consists of :

- a collarette
- a cylindrical body
- a rounded end

This plug may be introduced into the punctum without prior dilation.

The PAINLESS PLUG / SNUG PLUG is packaged in individual sterile blister packs, 2 per box. It is sterilized with ethylene oxide. A green indicator shows that the product followed a sterilization cycle validated by FCI.

The PAINLESS PLUG / SNUG PLUG (A) is delivered in an extended position thereby elongating the plug to facilitate insertion. It is pre-loaded on a disposable PLUG INSERTER (B) fitted with a clip (C) which keeps the plug in the extended position.

Indications:

The PAINLESS PLUG / SNUG PLUG is used to occlude the punctum in the treatment of ocular dryness. The closure is more effective when both puncta of a single side are closed.

Mode of action:

The three parts of the PAINLESS PLUG / SNUG PLUG are as follows:

The collarette covers the punctal opening on the eyelid and prevents the plug from entering too deeply. The collarette is thin in order to avoid corneal irritation. The extended rounded end enters the vertical section of the canaliculus without requiring prior dilation of the punctum. When released from the PLUG INSERTER, the plug resumes its natural expanded shape which allows it to fit snugly in place and prevents plug pop-out. The cylindrical body that connects the collarette to the rounded end has a blocked channel to allow for the insertion tool.

Contraindications:

Insertion of the PAINLESS PLUG / SNUG PLUG is not advisable in the following cases:

- Impermeability of the lacrimal ducts
- Eyelid abnormalities with absence of nocturnal occlusion
- Senile conjunctival hyperlaxity of the inferior punctum
- Positive result to a Rose Bengal test beyond the area of the palpebral fissure
- Rupture of the punctum

Adverse effects:

As with any type of surgery, there are risks linked to the insertion procedure of the plug and/or the evolution of the initial pathology. Possible complications associated with the insertion of a PAINLESS PLUG / SNUG PLUG include, but are not limited to the following:

Intraoperative complications (during the insertion of a PAINLESS PLUG / SNUG PLUG):

- Pain due to the dilation of the punctum
- Rupture of the punctum

Postoperative complications (once the PAINLESS PLUG / SNUG PLUG is in place):

- Temporary puritus that disappears after about a week
- Epiphora
- Chronic irritation
- Expulsion of the PAINLESS PLUG / SNUG PLUG
- Canalicular stenosis that may arise after the loss of a PAINLESS PLUG / SNUG PLUG
- Pyogenic granuloma

Unexpected side effects and complications that may be attributed to PAINLESS PLUGS / SNUG PLUGS must be reported to FCI.

Advice for use :

Dilation: Prior dilation of the punctum is not necessary. However, if done, any dilatation should be carried out with the dilating end of the applicator. Local anesthesia is recommended.

Insertion:

Insert the PAINLESS PLUG / SNUG PLUG with the PLUG INSERTER (fig. 1)

Once the PAINLESS PLUG / SNUG PLUG is in place in the punctum, (fig. 2),

Press on the handles (A) of the PLUG INSERTER until you feel a "click" which retracts the insertion needle and frees the PAINLESS PLUG / SNUG PLUG from the clip.

You must press the clip fully until the stop for complete extraction of the needle.

Removal: Removal is painless. Use forceps and pull on the cylindrical body just below the collarette.

Warnings and precautions for use:

The PAINLESS PLUG / SNUG PLUG must be removed from its blister pack and handled in aseptic conditions. Before use, the individual pack preserving the product's sterility should be checked to make sure it is intact. PAINLESS PLUGS / SNUG PLUGS are single-use products and must not be resterilized. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety. Do not use the product if the indicator is not green as it might endanger the patient. They should be stored at room temperature and should not be used after the expiration date shown on the package.

The guaranteed maximum residual rate of ethylene oxide the device can release to the patient is 4 µg. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

An implant card is available in the product's packaging. It is intended to be provided to the patient and allows traceability of the implanted product.

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