REF. 2138 307 86 C

FR	NOTICE D'UTILISATION - MONOKA de FAYET (Mandrin de Crawford)	2
EN	INSTRUCTION FOR USE – MONO-CRAWFORD / MONOKA of FAYET (Guide of Crawford)	4
DE	GEBRAUCHSANWEISUNG - MONOKA NACH FAYET (ZUGANG NACH CRAWFORD)	6
IT	ISTRUZIONI PER L'USO - MONOKA di FAYET (Mandrino di Crawford)	8
ES	INSTRUCCIONES DE USO - SONDA MONOKA de FAYET (Mandrill de Crawford)	10
PT	INSTRUÇÕES DE USO - MONOKA de FAYET (Mandrill de Crawford)	12
NL	GEBRUIKSAANWIJZING - MONOKA van FAYET (Mandrijn van Crawford)	14
DA	BRUGSVEJLEDNING - FAYETS MONOKA (Crawfords indføringsinstrument)	16
FI	KÄYTTÖOHJEET - FAYETIN MONOKA (CRAWFORD-OHJAIN)	18
SE	BRUKSANVISNING - MONOKA AV FAYET (LEDARE AV CRAWFORD)	19
TR	KULLANIM - FAYET MONOKA (Crawford Kılavuzu)	20
ZH	使用说明 - FAYET 的 MONOKA (Crawford 导引器)	22
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ – MONO-CRAWFORD / MONOKA of FAYET (Οδηγός Crawford)	23
PL	INSTRUKCJA UŻYTKOWANIA – STENTY MONO-CRAWFORD / MONOKA firmy FAYET (Prowadnica Crawford)	25
SL	NAVODILA ZA UPORABO – MONO-CRAWFORD/MONOKA of FAYET (CRAWFORDOVO VODILO)	27
SK	NÁVOD NA POUŽITIE – MONO-CRAWFORD / MONOKA of FAYET (vodidlo na Crawfordovu intubáciu)	29

Figures 31

Symbols 32

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MONO-CRAWFORD / MONOKA of FAYET (Guide of Crawford)

Description and packaging :

The MONO-CRAWFORD / MONOKA of FAYET (GUIDE OF CRAWFORD) is a monocalicular naso-lacrimal intubation device that is anchored at the punctum by plug. It is composed of :

a silicone tube (A)

a head (B) (plug) that ensures anchoring of the stent. This head is composed of the following three parts :

- an oval collarette
- a hollow body linking the collarette to the bulb
- a bulb attached to the silicone tube

a flexible steel guide (C)

The MONO-CRAWFORD / MONOKA OF FAYET (GUIDE OF CRAWFORD) is available in several models with different collarette sizes.

The MONO-CRAWFORD / MONOKA OF FAYET (GUIDE OF CRAWFORD) is also available with a polyvinylpyrrolidone (PVP) coating which improves its wettability and may enhance tear drainage while implant is in place.

The MONO-CRAWFORD / MONOKA OF FAYET (GUIDE OF CRAWFORD) is supplied sterile. A green indicator shows that the product followed a sterilization cycle validated by FCI. The MONO-CRAWFORD / MONOKA OF FAYET (GUIDE OF CRAWFORD) is sterilized using ethylene oxide in its final double packaging.

Indications :

Monocalicular nasolacrimal stents (MONO-CRAWFORD / MONOKA OF FAYET (GUIDE OF CRAWFORD)) are indicated in initial surgery in the treatment of monocalicular lacerations and stenoses requiring intubation. They are particularly useful in:

- congenital lacrimal duct obstructions
- canalicular lacerations
- after treatment of canaliculitis with concretions, when curettage of a canaliculus has been performed
- after dacryocystorhinostomy (conventional or laser)

Mode of action :

The MONO-CRAWFORD / MONOKA OF FAYET (GUIDE OF CRAWFORD) stents are used in naso-lacrimal duct intubation. The silicone tube takes the shape of the targeted region. In case of canalicular lacerations the silicone tube guides wound healing and reduces post-traumatic shrinkage. The plug-like shape of the head of the stent keeps the stent from migrating and from expelling.

Contraindications :

The MONO-CRAWFORD / MONOKA OF FAYET (GUIDE OF CRAWFORD) probes are contraindicated in the following cases:

- bicanalicular pathologies
- injuries of the lacrimal meatus
- need for prolonged intubation

Adverse side effects :

As in any type of surgery, there are risks linked to the material or to the evolution of the initial pathology. Potential complications associated with the implantation of a MONO-CRAWFORD / MONOKA OF FAYET (GUIDE OF CRAWFORD) stent include, but are not limited to the following:

Intraoperative complications :

- false passages
- separation of the silicone tube from the steel guide

Postoperative complications :

- | | |
|----------------------------------|--|
| - conjunctival or nasal pruritus | - Canalicular stenosis requiring stricturotomy |
| - nasal or caruncular irritation | - synechia of the nasal mucosa |
| - induced mucocele | - canaliculitis |
| - stent expulsion or loss | - keratitis or ulceration by rubbing of the collarette |
| - stent migration | |

Unexpected side effects and complications related to the MONO-CRAWFORD / MONOKA OF FAYET (GUIDE OF CRAWFORD) must be reported to FCI.

Advice for use :

The junction of the silicone to metal is a fragile area. It is recommended that the steel guide be held below this area while handling the stent.

Intubation technique:

The MONO-CRAWFORD / MONOKA OF FAYET (GUIDE OF CRAWFORD) requires the use of the Crawford Hook (S1.1275) which is specially designed to facilitate retrieval of the steel guide from the nose.

Advance the steel guide from the superior canaliculus into the floor of the nasal fossa.

Use a Crawford hook to establish metal-on-metal contact with the steel guide and to catch the olive tip.

Gently pull the hook with the steel guide out of the nose allowing the silicone stent to penetrate into the lacrimal ducts.

Use PLUG INSERTER (S1.3090) to seat the trailing plug in the punctum.

Warnings and precautions for use :

The MONO-CRAWFORD / MONOKA OF FAYET (GUIDE OF CRAWFORD) stents must be removed from their packaging and handled in aseptic conditions. Before use, the individual pack preserving the product's sterility should be checked to make sure it is intact. The MONO-CRAWFORD / MONOKA OF FAYET (GUIDE OF CRAWFORD) stents are single-use products and must not be resterilized. Reuse and/or reesterilization may compromise device performance, which could cause serious harm to the patient's health and safety. Do not use the product if the indicator is not green as it might endanger the patient. They should be stored at room temperature and must not be used after the expiration date shown on the package.

The guaranteed maximum residual rate of ethylene oxide at the time the device was placed on the market was 5.3 µg per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

An implant card is available in the product's packaging. It is intended to be provided to the patient and allows traceability of the implanted product.

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