REF. 2138 308 59 E

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SELF-RETAINING / AUTOSTABLE BICANALICULUS INTUBATION SET II WITH DILATOR

Description :

Each SELF-RETAINING / AUTOSTABLE BICANALICULUS INTUBATION SET II is composed of two silicone anchors connected to a silicone body delivered pre-mounted on two guides and packaged with a disposable dilator (Fig.1)

The SELF-RETAINING / AUTOSTABLE BICANALICULUS INTUBATION SET II is available in three different lengths :

25 mm	model : S1.1390
30 mm	model : S1.1391
35 mm	model : S1.1392

The sizer S1-1289 is used to determine the appropriate length (Fig.2) :

- Where relevant, treat the stenosis beforehand, and choose the appropriate length as follows :
 - Introduce the sizer into the inferior canaliculus,
 - Advance up to the bony contact while pulling gently on the eyelid,
 - Note location of the punctal opening on the sizer

Packaging :

SELF-RETAINING / AUTOSTABLE BICANALICULUS INTUBATION SET II is supplied sterile. It is sterilized by ethylene oxide. A green indicator shows that the product followed a sterilization cycle validated by FCI. Do not use the product if the indicator is not green. SELF-RETAINING / AUTOSTABLE BICANALICULUS INTUBATION SET II is sterilized in its final double-packaging to make it easier to handle in aseptic conditions.

Indications :

Bicanalicular intubation is indicated in treatments of epiphora in adults. The indications for bicanalicular intubation performed with the SELF-RETAINING / AUTOSTABLE BICANALICULUS INTUBATION SET II are :

- Punctal stenosis
- Canalicular stenosis within the lacrimal drainage system

Mode of action :

Bicanalicular intubations are used for stenting the canaliculi.
The silicone tube acts as a conformer and enables drainage of tears by capillarity.

Warning and contraindications :

The use of SELF-RETAINING / AUTOSTABLE BICANALICULUS INTUBATION SET II is contra-indicated in the following cases :

- Epiphora in infants
- Dacryolith
- Tumor of the lacrimal sac
- Other causes of epiphora that do not have a lacrimal origin, including, but not limited to, hypersecretion of tears (e.g. dry eyes, blepharitis, allergic conjunctivitis, corneal ulceration, etc.); failure of the lacrimal pump (e.g. facial paralysis)

The SELF-RETAINING / AUTOSTABLE BICANALICULUS INTUBATION SET II has not been studied for other conditions, such as previous eyelid and/or lacrimal surgery, a lump overlying or involving the punctum and/or other parts of the tear drainage system, long complete upper lacrimal system obstruction (canaliculi and common canaliculus) on diagnostic probing, or nasolacrimal duct stenosis or obstruction on irrigation testing.

Adverse effects :

Potential complications associated with the implantation of a SELF-RETAINING / AUTOSTABLE BICANALICULUS INTUBATION SET II include, but are not limited to the following ones :

- conjunctival pruritus
- exteriorization or loss of the silicone probe
- reversible shrinkage of the palpebral fissure
- Migration of the tube into the lacrimal passages if breakage of the device occurs
- caruncular irritation
- canaliculitis
- induced mucocele

The application of eye drops containing methylene blue can cause staining of the tube. This has no effect on the patient's tolerability of the stent. Unexpected side effects and complications related to the SELF-RETAINING / AUTOSTABLE BICANALICULUS INTUBATION SET II must be reported to FCI.

Recommendations for use :

See preparation and insertion below.

It is advisable to point out to the patient the interpalpebral loop following the procedure.

Preparation :

- Dilate each punctum using the disposable dilator included with the SELF-RETAINING / AUTOSTABLE BICANALICULUS INTUBATION SET II.
- Probe the inferior canaliculus and the superior canaliculus with a standard lacrimal probe to open any blockages.
- Pre-dilate the inferior and superior puncta using the disposable dilator included with the SELF-RETAINING / AUTOSTABLE BICANALICULUS INTUBATION SET II,
- Measure the canalicular length using the disposable sizer by introducing the sizer into the inferior canaliculus and advancing the sizer up to the bony contact, while pulling gently on the eyelid. Note the location of the punctal opening on the sizer and select the appropriate size stent (5,6 ou 7), as shown in Figure 2.
- If necessary, repeat dilation using the disposable dilator packaged with the SELF-RETAINING / AUTOSTABLE BICANALICULUS INTUBATION SET II.
- Insert the SELF-RETAINING / AUTOSTABLE BICANALICULUS INTUBATION SET II, following the instructions below.

Insertion :

- Hold the stent by the guides.
- With the help of the guide, place one anchor in the upper punctum while gently pulling on the eyelid
- Instruct the patient to look downward and carefully insert the stent along the axis of the canaliculus to a point beyond the central marker on the silicone body ; then remove the guide by immobilizing the stent body with a smooth-tipped forceps, such as monofilament tying forceps.
- Instruct the patient to look upward and insert the second anchor into the lower duct until the central marker is centrally positioned between the two puncta ; then remove the guide by immobilizing the silicone tube with the forceps.
- Have the patient close both eyes to bring the two puncta together, and gently massage the lacrimal sac to help position the stent properly.
- Verify the proper positioning of the stent by asking the patient to blink several times.

See Figure 3 below, for a properly inserted SELF-RETAINING / AUTOSTABLE BICANALICULUS INTUBATION SET II that has been advanced through the

canaliculi until the winglets have been deployed (A) and the central mark positioned mid-way between the puncta (B).

Removal :

- Grasp, with firm but gentle pressure, the visible loop with smooth-tipped forceps.
- With firm, gentle pressure, remove one end of the stent by pulling out from the punctum.
- Grasp the second end of stent near the other punctum and firmly remove the entire stent.

Precautions for use :

SELF-RETAINING / AUTOSTABLE BICANALICULUS INTUBATION SET II must be removed from its packaging and handled in aseptic conditions. Before use, the individual pack preserving the product's sterility should be checked to make sure it is intact. SELF-RETAINING / AUTOSTABLE BICANALICULUS INTUBATION SET II is a single-use product and must not be re-sterilized. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety. The stent should be stored between 41°F to 113°F and must not be used after the expiration date shown on the package.

The guaranteed maximum residual rate of ethylene oxide at the time the device was placed on the market was 5.3 µg per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

An implant card is available in the product's packaging. It is intended to be provided to the patient and allows traceability of the implanted product.

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