REF. 2138 309 01 D

FR	NOTICE D'UTILISATION - IBCN DE RITLENG®	2
EN	INSTRUCTION FOR USE - RITLENG® LACRIMAL INTUBATION SET	3
DE	GEBRAUCHSANWEISUNG - RITLENG® SCHLAUCHSET FÜR DIE TRÄNENWEGSINTUBATION	4
IT	ISTRUZIONI PER L'USO - KIT DI INTUBAZIONE LACRIMALE RITLENG®	5
ES	INSTRUCCIONES DE USO - IBCN DE RITLENG®	6
PT	INSTRUÇÕES DE USO - IBCN DE RITLENG®	7
NL	GEbruIKSAANWIJZING - IBCN VAN RITLENG®	8
DA	BRUGSVEJLEDNING – RITLENG® LACRIMALINTUBATIONSSÆT	9
FI	KÄYTTÖOHJEET - RITLENG®-LETKUSETTI KYYNELTIEHYEEN SONDEERAUKSEEN	10
SE	BRUKSANVISNING - RITLENG® SLANGSET FÖR INTUBATION AV TÄRKANA-LERNA	11
TR	KULLANIM - RITLENG® LAKRIMAL ENTÜBASYON SETİ	12
ZH	使用说明 - 泪道引流管 (Ritleng 型)	13
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - ΣΕΤ ΔΑΚΡΥΪΚΗΣ ΔΙΑΣΩΛΗΝΩΣΗΣ RITLENG®	15
PL	INSTRUKCJA UŻYTKOWANIA – ZESTAW DO INTUBACJI KANAŁÓW ŁZOWYCH RITLENG®	16
SL	NAVODILA ZA UPORABO – KOMPLET ZA LAKRIMALNO INTUBACIJO RITLENG®	17
SK	NÁVOD NA POUŽITIE - SÚPRAVA NA LAKRIMÁLNU INTUBÁCIU RITLENG®	18

Figure.....	19
Symbols.....	20

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EN

Instructions for use

Date of revision of the instructions for use : 04/2016

RITLENG® LACRIMAL INTUBATION SET**Description :**

The RITLENG® LACRIMAL INTUBATION SET is a nasal bicanalicular intubation for the lacrimal ducts. It consists of a silicone tube (external diameter: 0.64 mm) which is connected at each extremity with a PEEK thread guide. Each PEEK thread guide has a wide-diameter portion (0.39 mm) and a narrow-diameter portion (0.22 mm). The silicone tube may be coated with polyvinylpyrrolidone (PVP) to improve its wettability (PVP RITLENG® LACRIMAL INTUBATION SET).

Packaging :

The RITLENG® LACRIMAL INTUBATION SET is supplied sterile. It is sterilized by Ethylene Oxide. A green indicator shows that the product followed a sterilization cycle validated by FCI. Do not use the product if the indicator is not green. The RITLENG® LACRIMAL INTUBATION SETS are sterilized in their final double-pouch packaging to make them easier to handle in aseptic conditions.

Indications :

Nasal bicanalicular intubation is indicated in the treatment of epiphora in infants and adults:

- Canalicular pathologies (stenoses, obstructions, lacerations),
- Imperforation of the nasolacrimal duct in the infant (Epiphora in the infant),
- Dacryocystorhinostomy (conventional or laser).

This intubation can be used in patients of all ages and in all types of lacrimal surgery where an intubation is indicated.

Mode of action :

The RITLENG® LACRIMAL INTUBATION SET is used to intubate the lacrimal ducts. The silicone tube acts as a conformer and enables drainage of tears by capillarity. In cases of canalicular lacerations, the silicone tube guides wound healing and reduces post-traumatic shrinkage. The advantages of the RITLENG® LACRIMAL INTUBATION SET over conventional bicanalicular probes are the following :

Removing a PEEK thread is less traumatic to the nasal mucosa, less haemorrhagic and less painful under local anesthesia, than removing a metal guide. Locating a PEEK buckle that lies wide and flat in the nose is much easier than locating a metal guide that comes out in the inferior meatus, often difficult to reach.

Contraindications :

As with all other forms of intubation, a visual monitoring of the nasal cavity is recommended with the help of intranasal illumination or a nasal endoscope (0° or 30°) which makes it easier to locate the thread guide, particularly if it takes a posterior direction. It is recommended not to remove the PEEK when the guide probe is in the nasolacrimal duct as it can get hooked in the groove of the probe becomes fragile and tears.

There is no contra-indication specific to this intubation.

Side effects :

As in any type of surgery, there are risks linked to the material or to developments of the initial pathology. Potential complications associated with the implantation of a RITLENG® LACRIMAL INTUBATION SET are similar to those associated with conventional nasal bicanalicular intubations, the most frequent being :

Intraoperative complications :

- false passages
- separation of the tube and the PEEK thread
- difficulty in locating the thread in the nose (hence the importance of endoscopy).

Postoperative complications (rare):

- conjunctival pruritus
- irritation in the internal angle of the eye
- sneezing
- stricturotomy
- epistaxis (unusual)
- canaliculitis
- exteriorization of the probe

Unexpected side effects and complications related to the RITLENG® LACRIMAL INTUBATION SET must be reported to FCI.

Instructions for use :

Using the RITLENG® LACRIMAL INTUBATION SET requires use of the RITLENG® PROBE (S1-1460) which is specially designed for this purpose. The intubation technique presents the following particularities:

Probing of the first canaliculus towards the nasolacrimal duct, up to the floor of the nasal fossa.

After this probing with the RITLENG® PROBE and guiding its plateau to direct the inferior opening towards the front of the nose, the PEEK thread guide is introduced in the RITLENG® PROBE so as to obtain a wide loop which will lie flat on the nasal fossa.

It is strongly recommended that a nasal endoscope and fine endonasal forceps be used for easy recovery of the PEEK thread.

In the absence of endoscopy, take a small foam hook, seek contact with the inferior extremity of the metal probe and rotate the guide probe towards the rear.

Then gently retract the probe and as soon as contact with the metal is lost, pull gently on the hook and remove the thread.

The PEEK thread (wide section) is removed from the nose and the metal probe is extracted from the lacrimal ducts.

The RITLENG® PROBE is then removed from the intubation by pending the narrow section (0.22 mm) of the PEEK thread guide out through the slit.

The same procedure is followed to probe the second canaliculus to achieve a nasal bicanalicular intubation, with an internal angle loop that is neither too loose nor too tight.

It is recommended that the patient be shown the interpalpebral loop after the procedure.

Precautions for use :

RITLENG® LACRIMAL INTUBATION SETS must be removed from their pouch and handled in aseptic conditions. Before use, the individual pack preserving the product's sterility should be checked to make sure it is intact. RITLENG® LACRIMAL INTUBATION SETS are single-use products and must not be resterilized. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety. They should be stored at room temperature and must not be used after the expiration date shown on the package.

The guaranteed maximum residual rate of ethylene oxide at the time the device was placed on the market was 5.3 µg per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

An implant card is available in the product's packaging. It is intended to be provided to the patient and allows traceability of the implanted product.

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