REF. 2138 309 02 E

FR	NOTICE D'UTILISATION - MONOKA ®	2
EN	INSTRUCTION FOR USE - MONOKA ®	4
DE	GEBRAUCHSANWEISUNG - MONOKA ®	6
IT	ISTRUZIONI PER L'USO - MONOKA ®	8
ES	INSTRUCCIONES DE USO - MONOKA ®	10
PT	INSTRUÇÕES DE USO - MONOKA ®	12
NL	GEBRUIKSAANWIJZING - MONOKA ®	14
DA	BRUGSVEJLEDNING - MONOKA ®	16
FI	KÄYTTÖOHJEET - MONOKA ®	18
SE	BRUKSANVISNING - MONOKA ®	20
TR	KULLANIM - MONOKA ®	22
ZH	使用说明 - 泪道引流管 (MONOKA 单管型)	24
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - MONOKA ®	27
PL	INSTRUKCJA UŻYTKOWANIA – MONOKA ®	29
SL	NAVODILA ZA UPORABO – MONOKA ®	31
SK	NÁVOD NA POUŽITIE - MONOKA ®	33

Symbols 35

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EN

Instructions for use
Date of revision of the instructions for use : 06/2016

MONOKA®

Description :

The MONOKA® stents are monocanalicular-nasal or monocanalicular intubation sets with stable meatic fixation, that are anchored at the punctum by plug. They are composed of a silicone tube with the following at each of its extremities:

- a head (plug) that ensures "stable meatic fixation" of the probe. This head is composed of the following three parts :
 - an oval collarette
 - a hollow body linking the collarette to the bulb,
 - a bulb fixed to the silicone tube,
- and depending on the type of MONOKA® :
 - either a malleable steel guide (FAYET-BERNARD MONOKA®),
 - or a PEEK thread guide (FAYET-BERNARD MONOKA® RITLENG® TECHNIQUE).
- No guide in case of the MINI-MONOKA® .

There are different models of the MONOKA® stents which differ in collarette or steel guide sizes. The MONOKA® stents may be coated with polyvinylpyrrolidone (PVP) to improve their wettability.

Packaging :

The MONOKA® stents are supplied sterile. They are sterilized by ethylene oxide in their final packaging. A green indicator shows that the product followed a sterilization cycle validated by FCI. Do not use the product if the indicator is not green. The MONOKA® are sterilized in their final double-blister or double-pouch packaging to make them easier to handle in aseptic conditions.

Indications :

The nasolacrimal monocanalicular stents (MONOKA®) are indicated in initial surgery process in the treatment of monocanalicular laceration and stenoses requiring intubation. They are particularly used in:

- congenital nasolacrimal duct obstructions
- canalicular laceration
- after treatment of canaliculitis with concretions, when curettage of a canaliculus has been performed.
- after dacryocystorhinostomy (conventional or laser)

True monocanalicular stents (MINI-MONOKA®) are mainly indicated in traumatology with injuries of the outer 2/3 of one of the canaliculi. They can also be indicated in stenoses of the inferior canaliculus (post-infectious or not).

Mode of action :

The MONOKA® stents are used in nasal monocanalicular or in monocanalicular intubation. The silicone tube takes the shape of the targeted region. In case of canalicular lacerations the silicone tube guides wound healing and reduces post-traumatic shrinkage. The shape of the head of the stent (plug) keeps the probe from migrating and from expelling.

Warning and contraindications :

The MONOKA® stents are contra-indicated in the following cases:

- bicanalicular pathologies.
- injuries of the lacrimal meatus.
- need for prolonged intubation.

The silicone-metal or silicone-PEEK junction is a fragile area. It is recommended that the steel guide or the thread guide be held in front of this area while handling the intubation.

Adverse side effects :

As in any type of surgery, there are risks linked to the material or to developments of the initial pathology. Potential complications associated with the implantation of a MONOKA® stent include, but are not limited to the following :

Intraoperative complications :

- false passages
- separation of the silicone tube and the steel guide or of the silicone tube and the PEEK thread guide.

Postoperative complications :

- conjunctival or nasal pruritus
- exteriorization or loss of the probe
- synechia of the nasal mucosa
- nasal or caruncular irritation
- migration of the probe
- canaliculitis
- induced mucocele
- stricturotomy
- keratitis or ulceration by rubbing of the collarette

The application of eye drops containing methylene blue can cause staining of the tube. This has no effect on the patient's tolerability of the probe.

Unexpected side effects and complications related to the MONOKA® must be reported to FCI.

Instructions :

It is recommended, when appropriate, to perform intubation by the superior duct.

Intubation technique :

Introduce the silicone tube into the lacrimal duct until the plug is at the punctal opening. Use a PLUG INSERTER 0.3 (S1.3090 or S1.3090U) to place the plug.

Particular Cases :

Using the RITLENG® TECHNIQUE MONOKA requires use of the RITLENG® PROBE (S1.1460 OR S1.1460U) which is specially designed for this purpose.

The intubation technique requires the following steps:

Prepare the nasal fossa with a vasoconstrictive packing

Dilate punctum and perform a preliminary probing of lacrimal ducts with a Bowman probe

Introduce the RITLENG® PROBE into the canaliculus and into the naso-lacrimal duct After contact with the nasal fossa floor, pull the probe out slightly (1 cm) to provide access for the PEEK thread to spread out in the nose.

Orient the probe disc to direct the inferior opening of the probe towards the nose entrance (push the disc back to direct the bottom of the probe towards the nostril)

Feed a generous supply of the PEEK thread into the probe to form a large buckle that will spread out in the nasal cavity

There are several ways to remove the PEEK thread from the nasal cavity :

Occasionally the thread ends up the entrance of the nasal fossa : it is easy to grasp it with forceps (a little aspiration may help) or with a blunt hook .

The use of a nasal endoscope is very helpful in removing the thread under visual control especially if the PEEK thread ends up in the rear of the nasal fossa, which sometimes happens. A fine endonasal forceps (FCI) facilitates the recovery of the thread.

When no endoscope is used, use the blunt RITLENG® HOOK to make the metal-to-metal contact with the probe at the junction of the lateral wall and the floor of the nasal fossa. When contact is made, turn the probe 180°, orienting it towards the rear (this rotation makes it easier to grasp the thread with the hook). Then slowly back the probe up the naso-lacrimal duct and just as the metal-to-metal contact is lost remove the hook from the nose pulling the PEEK thread

along with it.

When the wide portion of the PEEK thread starts to exit the nose, hold it there while removing the RITLENG® PROBE back out from the lacrimal ducts.

When the thinner portion of the thread enters the probe, separate it from the probe by sliding it out through the open slit (before the advance of remaining thicker portion of the thread and silastic tubing).

Use a PLUG INSERTER to seat the plug into the punctum (do not pull the plug into place).

Using the MINI-MONOKA® involves the following particularities:

Intubation of the external segment after careful dilation of the lacrimal meatus,

Identification of the internal segment of the canaliculus, after which the probe is introduced in this internal segment with the help of Paufrage forceps. The canaliculus is then sutured.

Precautions for use :

MONOKA® stents must be removed from their packaging and handled in aseptic conditions. Before use, the individual pack preserving the product's sterility should be checked to make sure it is intact. MONOKA® PROBES are single-use products and must not be resterilized. Reuse and/or reesterilization may compromise device performance, which could cause serious harm to the patient's health and safety. They should be stored at room temperature and must not be used after the expiration date shown on the package.

An implant card is available in the product's packaging. It is intended to be provided to the patient and allows traceability of the implanted product.

The guaranteed maximum residual rate of ethylene oxide at the time the device was placed on the market was 5.3 µg per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

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