REF. 2138 309 25 J

FR	NOTICE D'UTILISATION - MASTERKA®	2
EN	INSTRUCTION FOR USE - MASTERKA®	4
DE	GEBRAUCHSANWEISUNG - MASTERKA®	6
IT	ISTRUZIONI PER L'USO - MASTERKA®	8
ES	INSTRUCCIONES DE USO - MASTERKA®	10
PT	INSTRUÇÕES DE USO - MASTERKA®	12
NL	GEBRUIKSAANWIJZING - MASTERKA®	14
DA	BRUGSVEJLEDNING - MASTERKA®	16
FI	KÄYTTÖOHJEET - MASTERKA®	18
SE	BRUKSANVISNING - MASTERKA®	20
TR	KULLANIM - MASTERKA®	22
ZH	使用说明 - MASTERKA®	24
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - MASTERKA®	26
SL	INSTRUKCJA UŻYTKOWANIA – MASTERKA®	28
PL	NAVODILA ZA UPORABO – MASTERKA®	30
SK	NÁVOD NA POUŽITIE - MASTERKA®	32

Figures **34**

Symbols **35**

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MASTERKA®

Description and packaging:

MASTERKA® is a device with self-retaining punctal fixation. It consists of a single unit including a silicone tube with a fixation head and an introducer (see diagram). There are 3 models.

- The silicone tube is cylindrical with one end closed and the other end closed by a fixation head.
- The fixation head consists of three parts:
 1. A bulb intended to prevent the spontaneous extrusion of the tube
 2. A neck which will be in contact with the punctum.
 3. A collarette that will lay flush on the eyelid margin and prevent migration of the fixation head. The fixation head is molded with an opening to make room for a plug inserter.
- The introducer facilitates the insertion of the MASTERKA®. The longer portion of the introducer is placed through the opening and into the lumen of the silicone tube. It is completely removed once the intubation of the lacrimal passages has been completed.

The free end of the introducer is covered with a cylindrical sleeve to serve as a handle and facilitate manipulation.

MASTERKA® is supplied sterile and pre-mounted onto its introducer. It is sterilized by ethylene oxide. A green indicator shows that the product followed a sterilization cycle validated by FCI. The MASTERKA® has been sterilized in its final double packaging to facilitate handling in a sterile environment.

Indications:

The MASTERKA® is a silicone intubation device indicated for use in the treatment of congenital lacrimal duct obstructions (stenosis of the valve of Hasner) in patients 12 months and older.

Mode of action:

The silicone stent acts as a conformer. The plug-like shape of the head of the stent keeps the stent from migrating and from expulsing.

Contra-indications :

- a) The application of eye drops containing methylene blue can cause staining of the tube. This has no effect on the patient's tolerance of the stent.
- b) The use of the MASTERKA® is not recommended and/or is contraindicated in the following cases:
- Bicanalicular pathologies
 - Stenosis of the common canaliculus
 - Congenital or acquired enlargement of the lacrimal meatus (stricturotomy ...)
 - Tumor of the lacrimal sac
 - Adult nasolacrimal stenosis
 - Dacryolith
 - Complex congenital nasolacrimal obstruction

Adverse side effects:

Just as in any type of surgery there exist risks associated with the material and/or the evolution of the initial pathology. The potential complications associated with the insertion of the MASTERKA® tube include :

Complications which may arise during the application of the tube:

- Rupture of the lacrimal meatus due to excessive dilation.
- False passages or perforation of the tube with the introducer, which may result in canalicular trauma or bleeding.

Post-operative complications:

- Conjunctival caruncular pruritis or half-moon folding
- Keratitis or ulceration caused by the rubbing of the ring
- Extrusion or loss of the tube
- Migration of the tube into the lacrimal passages
- Nasal Pruritis
- Canaliculitis
- Resulting mucocele
- Synechiae of the nasal mucous membranes

MASTERKA® tubes may cause Secondary effects and complications that must be reported to FCI.

Instructions for use:

Exploration:

Insertion of the MASTERKA® must be preceded by exploration of the lacrimal passages with a surgical instrument such as a Bowman probe. The result of this exploration will determine which stent to use and which model in the case of a Masterka. The location and severity of obstruction is identified and if the obstruction is at the valve of Hasner, the Masterka is appropriate for use.

During the exploration, it is important to avoid false passages and to check for proper placement by carefully introducing a longer metallic probe with a blunt tip under lower nasal meatus in search of the metal-to-metal contact of the 2 probes (Fig.1)

Selection of MASTERKA® model:

In order to determine the appropriate Masterka model, the distance between the punctum and the nasal floor is measured. The MASTERKA® should be equal to or shorter than this distance, but slightly longer than the distance between the punctum and the obstruction (Fig.2).

FCI's measuring device S1.1289 is used to select the appropriate Masterka length.

Insert the measuring device until it reaches the nasal cavity floor. Confirm proper placement with another blunt-tipped metallic probe by achieving metal-to-metal contact with the measuring device (Fig.1)

Locate the marking on the measuring device that corresponds to the punctal opening, (Fig.1)

Select the correct MASTERKA® model according to the markings on the measuring device (1,2 or 3) of Fig.3.

Stent insertion:

- Once the necessary conditions have been determined (proper therapeutic need, absence of contraindications, anaesthesia based on patient age and condition), the lacrimal duct is dilated.
- The MASTERKA®stent is inserted into the canaliculus which has been gently dilated.
- The MASTERKA®stent is inserted until it makes contact with the bone; the introducer is then rotated in order to catheterise the lacrimal sac and the nasolacrimal passage. This movement is the same as that used in the exploration of the lacrimal passages. This also positions the fixation head in contact with the punctal opening (Fig. 4).
- To avoid false passages and to check proper placement of the MASTERKA®, to introduce carefully in the same way the probe than used to detect false passages during the exploration, in search of the contact between the probe and stent.
- Removal of the introducer: The fixation head should be firmly held against the punctal opening while at the same time the introducer is fully extracted from the silicone tube (Fig. 4).
- The fixation head is then seated into the vertical canaliculus with the help of a plug inserter FCI (S1.3090).
- The collarette should be firmly placed flush against the eyelid margin. The fixation head neck should no longer be visible.

Warnings and precautions for use :

ATTENTION: TECHNICALS ERRORS :**False passages (Fig.1, Fig.6):**

- the detection of false passages during the instrumental exploration step and during the use of the measuring device is realized by searching metal-to-metal contact between exploratory probe or the measuring device and a second blunt-tipped probe under lower nasal meatus.
- the detection of false passages after the insertion of the MASTERKA® is realized by searching metal-to-MASTERKA® contact between Masterka® and the second blunt-tipped probe under lower nasal meatus.

Stent too long (Fig.7):

- the end of the stent rubs against the floor of the nasal cavity but the fixation head remains above the punctum. In this case, select a shorter MASTERKA®

Wrong position (Fig.8):

- The introducer is fully removed and the fixation head remains 2 or 3 mm above the punctum. One should not try to force the stent in further with the applicator, instead check the length of the stent and restart the procedure from the beginning.

Stent perforation (Fig.9, Fig.10):

during the insertion of the stent MASTERKA® and in the case of complex stenosis, it is possible to perforate the silicone tube with the introducer thus preventing the tube from passing through the obstruction (Fig.9). This is confirmed if the distance between the fixation head and the introducer's handle is too short (Fig.10) and if upon removal of the introducer, the stent backs out along with it.

The MASTERKA® stent must be removed from its packaging and handled under sterile conditions. Check that the individual sterile protection is intact before use. The MASTERKA® tube is a single-use system and should not be re-sterilized. Reuse and/or re-sterilization may compromise device performance, which could cause serious harm to the patient's health and safety. Do not use the product if the indicator is not green as it might endanger the patient. The MASTERKA® tube should be stored at room temperature and should not be used after the use-by date indicated on the package.

The guaranteed maximum residual rate of ethylene oxide at the time the device was placed on the market was 5.3 µg per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

An implant card is available in the product's packaging. It is intended to be provided to the patient and allows traceability of the implanted product.

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