**REF. 2138 309 34 B**

FR	NOTICE D'UTILISATION - LENTILLES À USAGE UNIQUE EN PMMA POUR VITRECTOMIE	2
EN	INSTRUCTION FOR USE - DISPOSABLE PMMA VITRECTOMY LENS	3
DE	GEBRAUCHSANWEISUNG - VITREKTOMIELINSE AUS PMMA FÜR DEN EINMALGEBRAUCH	4
IT	ISTRUZIONI PER L'USO - LENTI MONOUSO IN PMMA PER VITRECTOMIA	5
ES	INSTRUCCIONES DE USO - LENTE PMMA DESECHABLE PARA VITRECTOMÍA	6
PT	INSTRUÇÕES DE USO - LENTE DE VITRECTOMIA DESCARTAVEL EM PMMA	7
NL	GEBRUIKSAANWIJZING - VITRECTOMIELENS VAN PMMA VOOR EENMALIG GEBRUIK	8
DA	BRUGSVEJLEDNING - ENGANGS PMMA VITREKTOMILINSE	9
FI	KÄYTTÖOHJEET - KERTAKÄYTTÖINEN PMMA-VITREKTOMIALINSSI	10
SE	BRUKSANVISNING - VITREKTOMIELINSE AUS PMMA FÜR DEN EINMALGEBRAUCH	11
TR	KULLANIM - TEK KULLANIMLIK PMMA VITREKTOMİ LENSİ	12
ZH	使用说明 - 玻璃体切割术用一次性有机玻璃透镜	13
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - ΦΑΚΟΣ ΥΑΛΟΕΙΔΕΚΤΟΜΗΣ ΜΙΑΣ ΧΡΗΣΗΣ ΑΠΟ PMMA	14
PL	INSTRUKCJA UŻYTKOWANIA – JEDNORAZOWA SOCZEWKA DO WITREKTOMII PMMA	15
SL	NAVODILA ZA UPORABO – LEČA ZA VITREKTOMIJO PMMA ZA ENKRATNO UPORABO	16
SK	NÁVOD NA POUŽITIE - JEDNORAZOVÁ VITREKTOMICKÁ ŠOŠOVKA PMMA	17

**Symbols..... 18**

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## DISPOSABLE PMMA VITRECTOMY LENS

### Description :

The DISPOSABLE VITRECTOMY LENSES are made of polymethylmethacrylate (PMMA).

### Packaging :

The DISPOSABLE PMMA VITRECTOMY LENSES are supplied sterile. They are sterilized using ethylene oxide. A green indicator shows that the product followed a sterilization cycle validated by FCI. Do not use the product if the indicator is not green.


For ease of handling in aseptic conditions, the DISPOSABLE PMMA VITRECTOMY LENSES are presented in a double packaging.

### Indications :

The DISPOSABLE PMMA VITRECTOMY LENSES are used in the context of vitrectomy surgeries.

The lenses, used with a microscope, facilitate visualisation of the retina by cancelling the corneal dioptré.

A raised lip (with 3 stability pods for a PLANO-CONCAVE AUTOSTABLE LENS) facilitates the handling and positioning of the lens on the cornea using forceps. The top of this lip is left unpolished in order to distinguish it from the polished side that is placed on the cornea.

### Complications:

Any secondary effects or unexpected complications possibly linked to the use of the DISPOSABLE PMMA VITRECTOMY LENSES must be reported to FCI .

### Intructions for use :

To be used with a microscope.

The surfaces are polished and care must be taken to ensure that they are not damaged or dirtied.

Before applying the lens against the cornea, lubricate its concave side with a drop of ophthalmic gel.

### Warnings :

The lenses must be stored at room temperature.

The DISPOSABLE PMMA VITRECTOMY LENS should be removed from its packaging and handled in aseptic conditions. Before use, the individual pack preserving the product's sterility should be checked to make sure it is intact. The DISPOSABLE PMMA VITRECTOMY LENSES are single-use products and must not be resterilized. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety. It should be kept at room temperature and must not be used after the expiration date shown on the package.

The guaranteed maximum residual rate of ethylene oxide at the time the device was placed on the market was 5.3 µg per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

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