REF. 2138 310 96 C

FR	NOTICE D'UTILISATION - FCI NUNCHAKU ®	2
EN	INSTRUCTION FOR USE - FCI NUNCHAKU ®	4
DE	GEBRAUCHSANWEISUNG - FCI NUNCHAKU ®	6
IT	ISTRUZIONI PER L'USO - FCI NUNCHAKU ®	8
ES	INSTRUCCIONES DE USO - FCI NUNCHAKU ®	10
PT	INSTRUÇÕES DE USO - FCI NUNCHAKU ®	12
NL	GEBRUIKSAANWIJZING - FCI NUNCHAKU ®	14
DA	BRUGSVEJLEDNING - FCI NUNCHAKU ®	16
FI	KÄYTTÖOHJEET - FCI NUNCHAKU ®	18
SE	BRUKSANVISNING - FCI NUNCHAKU ®	20
TR	KULLANIM - FCI NUNCHAKU ®	22
ZH	使用说明 - FCI NUNCHAKU ®	24
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - FCI NUNCHAKU ®	26
PL	INSTRUKCJA UŻYTKOWANIA – FCI NUNCHAKU ®	28
SL	NAVODILA ZA UPORABO – FCI NUNCHAKU®	30
SK	NÁVOD NA POUŽITIE - FCI NUNCHAKU ®	32

Figures **34**

Symbols **35**

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EN

Instructions for use
Date of revision of the instructions for use: 12/2016

FCI NUNCHAKU®

Description and presentation:

The FCI NUNCHAKU® is a silicone self-retaining bicanaliculus intubation stent, mounted on metallic guide to facilitate insertion.

The stent ends are cone-shaped and one of the ends is colored to serve as a reference for the surgeon during implantation.

The middle part is flexible and the sides are rigid.

Both the stent and the guide have markings that act as references for the surgeon during implantation (refer to the attached intubation set diagrams).

The FCI NUNCHAKU® is supplied sterile and pre-mounted onto guides. It is sterilized using ethylene oxide. A green indicator shows that the product followed a sterilization cycle validated by FCI.

The FCI NUNCHAKU® is sterilized inside its final double-packaging to facilitate handling under aseptic conditions.

Indications:

The FCI NUNCHAKU® stents are indicated for use for nasolacrimal intubation in patients 12 months and older. Indications for nasolacrimal intubation performed with the FCI NUNCHAKU® are:

- canalicular pathologies (e.g., congenital or acquired stenosis, lacerations)
- during dacryocystorhinostomy
- congenital lacrimal duct obstruction

Mode of action:

The FCI NUNCHAKU® acts as a conformer, allowing tears to be drained by capillarity.

Contraindications:

The use of FCI NUNCHAKU® is contraindicated in the following cases:

- Epiphora in infants,
- Dacryolith,
- Tumor of the lacrimal sac,
- Other causes of epiphora that do not have a nasolacrimal origin, including, but not limited to, hypersecretion of tears (e.g. dry eyes, blepharitis, allergic conjunctivitis, corneal ulceration, etc.); or, failure of the lacrimal pump (e.g. facial paralysis).

Adverse side effects:

As in any type of surgery, there are potential risks linked to the material and/or to the progression of the initial pathology. Potential complications associated with the implantation of the FCI NUNCHAKU® intubation set include, but are not limited to, the following:

Intraoperative complications:

- False passages,
- Canalicular trauma.

Postoperative complications:

- | | |
|---|---|
| <ul style="list-style-type: none"> - Conjunctival or nasal pruritus, - Reversible shrinkage of the palpebral fissure, - Strictureotomy, - Canaliculitis, - Migration of the tube into the lacrimal passages if breakage of the intrapalpebral loop occurs, | <ul style="list-style-type: none"> - Nasal or caruncular irritation, - Exteriorization or loss of the silicone stent, - Synechia of the nasal mucosa, - Induced mucocele, - Infection, - Bleeding |
|---|---|

Unexpected side effects and complications related to the FCI NUNCHAKU® must be reported to FCI.

Recommendations for use:

Preparation:

Probe the lacrimal tear ducts with a standard lacrimal probe to open any blockages before implanting the FCI NUNCHAKU® stent.

Selecting the stent:

It is recommended to use a 90mm stent for children and a 105mm stent for adults.

Implanting the stent:

1. After probing the lacrimal ducts, insert the stent through the punctum and into the vertical part of the canaliculus, taking care to use only the guide to insert the stent. Grasping the stent along with the guide during insertion could result in perforation of the stent tip.
2. Guide the stent parallel to the eyelid margin. When the FCI NUNCHAKU® makes contact with the bone; the guide is then rotated in order to catheterise the lacrimal sac and the nasolacrimal passage. This movement is the same as that used in the exploration of the lacrimal passages. If a partial stenosis of the canaliculus or a major curvature of the nasolacrimal duct causes resistance to be felt while advancing the stent, remove the metallic guide from the stent a few millimeters to gain flexibility. Then, advance both the silicone stent and the metallic guide using delicate forceps. This will prevent the silicone stent from creating a false passage.
3. To check for proper placement of the NUNCHAKU®, a longer metallic probe with a blunt tip is carefully introduced under the lower nasal meatus in search of the metal-to-NUNCHAKU® contact.
4. Remove the guide by holding the silicone tube in place with delicate forceps and gently rotating the handle of the guide while pulling upwards until the guide is fully removed from the stent lumen.
5. The same procedure is followed to stent the second canaliculus in order to achieve a bicanalicular intubation. After both ends of the stent are inserted, make sure the interpalpebral loop connecting the stent is centered and neither too loose nor too tight.
6. Have the patient close both eyes to bring the two puncta together, and gently massage the lacrimal sac to help position the stent properly.
7. Check the proper positioning of the stent by asking the patient to blink several times and to verify that the tension on the intrapalpebral loop remains the same and the stent remains in the proper location.
8. It is recommended that the patient be shown the interpalpebral loop after the surgery, instructing the patient to avoid pulling on the loop, which could dislodge the stent.

Stent removal:

Removal of the silicone tubes is accomplished by locating the silicone tube between the upper and lower puncta, grasping firmly with tweezers or forceps, and pulling the tubing upward.

Warnings and precautions for use:

The application of eye drops containing methylene blue can cause staining of the tube. This has no effect on patient's tolerability of the stent.

The FCI NUNCHAKU® must be removed from its packaging and handled under aseptic conditions. Check the integrity of the individual pack before use.

The FCI NUNCHAKU® is a single use device that must not be re-sterilized. Re-use and/or re-sterilization could affect the device performances and compromise the patient's safety and health. Don't use the product if the sterilization indicator is not green, as it could endanger the patient.

The FCI NUNCHAKU® should be stored at room temperature and should not be used after the expiration date shown on the package.

Designer

Katsuaki KURIHASHI, M.D.

The guaranteed maximum residual rate of ethylene oxide at the time the device was placed on the market was 5.3 µg per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

An implant card is available in the product's packaging. It is intended to be provided to the patient and allows traceability of the implanted product.

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