FR	NOTICE D'UTILISATION – LACRIJET®	2
EN	INSTRUCTION FOR USE – LACRIJET®	4
DE	GEBRAUCHSANWEISUNG – LACRIJET®	6
IT	ISTRUZIONI PER L'USO – LACRIJET®	8
ES	INSTRUCCIONES DE USO – LACRIJET®	10
PT	INSTRUÇÕES DE USO – LACRIJET®	12
NL	GEBRUIKSAANWIJZING – LACRIJET®	15
DA	BRUGSVEJLEDNING – LACRIJET®	17
FI	KÄYTTÖOHJEET – LACRIJET®	19
SV	BRUKSANVISNING – LACRIJET®	21
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ – LACRIJET®	23
ET	KASUTUSJUHEND – LACRIJET®	25
LV	LIETOSANAS INSTRUKCIJA – LACRIJET®	27
LT	NAUDOJIMO INSTRUKCIJOS – LACRIJET®	29
PL	INSTRUKCJA UŻYTKOWANIA – LACRIJET®	31
SL	NAVODILA ZA UPORABO – LACRIJET®	34
SK	NAVOD NA POUZIVANIE – LACRIJET®	36
CS	NAVOD K POUZITI – LACRIJET®	38
RO	INSTRUCȚIUNI DE UTILIZARE – LACRIJET®	40
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА – LACRIJET®	42
HR	UPUTE ZA UPOTREBU – LACRIJET®	44
HU	HASZNALATI UTASITAS – LACRIJET®	46
NO	BRUKSANVISNING – LACRIJET®	48
TR	KULLANIM – LACRIJET®	50
AR	تعليمات الاستخدام – LACRIJET®	52
Illustrations .....		55
Symbols .....		56
<b>REF. 2138 312 23 E</b>		<b>DESIGNER: Bruno FAYET, M.D.</b>

FCI S.A.S.  
20-22, rue Louis Armand  
75015 PARIS – FRANCE

☎ : + 33 (0) 1 53 98 98 98

📠 : + 33 (0) 1 53 98 98 99

@ : fci@fci.fr

🌐 : www.fci.fr

Distributed in U.S.A. by:

FCI Ophthalmics Inc.  
30 Corporate Park Drive  
Suite # 310/320  
Pembroke, MA 02359

☎ : 800-932-4202

📠 : 781-826-9062

@ : info@fci-ophthalmics.com

🌐 : www.fci-ophthalmics.com

R<sub>x</sub> Only

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LACRIJET® is a mono-nasocanalicular stent of the Monoka® type, with self-retaining punctal fixation, pre-loaded into a handpiece to facilitate the introduction of the stent into place in the tear ducts. LACRIJET® comes in 7 models with stents from 15 to 50mm.

The silicone stent consists of a solid silicone tube with an external diameter of 0.64 mm, and a head for the self-retaining fixation of the stent. This head is made of the following three parts:

- a plug to prevent the stent from spontaneous expulsion (A),
- a hollow body that connects the collarette to the plug (B),
- a collarette that caps the punctal opening in order to prevent intra-canalicular migration of the stent (C).

A lumen passes through the fixation head to allow for a plug inserter (ref S1.3090) to be used.

The handpiece, onto which the stent is pre-loaded consists of 4 main parts:

- a body for gripping the device (D),
- a metal guide-tube on which the stent is loaded (E),
- a detachable unlocking part (F),
- a piston, inside the body that slides backwards so that the stent can be inserted (G).

The device is sterilized using irradiation in its final double blister packaging to make it easier to handle in aseptic conditions. A red indicator shows that the product followed a sterilization cycle validated by FCI. The device is for single use only and is not resterilizable.

#### Indications

LACRIJET® from 30 to 50mm is indicated in the treatment of congenital lacrimal duct obstructions (stenosis of the valve of Hasner) in patients 12 months and older.

LACRIJET® with small size (15 and 20mm) is indicated in the repair of canalicular lacerations.

This medical device is intended for use by healthcare professional only.

#### Contraindications

Utilization of LACRIJET® is contraindicated in the following cases:

- Bicanalicular pathologies
- Stenosis of the common canaliculus
- Tumor of the lacrimal sac
- Adult nasolacrimal stenosis
- Congenital or acquired enlargement of the lacrimal meatus (stricturotomy)
- Dacryolith
- Complex congenital nasolacrimal obstruction
- Filiforme nasal canal.

#### Warnings

The device is for single use and must not be resterilized. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety. Check packaging integrity before use, and do not use the device if package is damaged. Do not use the product if the indicator is not red, as it could endanger the patient. Do not use after the expiration date shown on the package.

The LACRIJET® implantable stent is made entirely of medical grade silicone. The device is usually left in place for 3 to 6 months.

Care should be given to patients with dacryocystitis and craniofacial abnormalities because of increased risks of complications.

An attention to eye-rubbing behavior must be given since this behavior might be linked with premature exteriorization and/or loss of the stent, with the associated consequences (complications and intubation failure). As a rule, a particular attention should be given to patients that may not follow directions (children and patient with neurologic disorders.)

In case the user identifies faulty labelling information, it shall be relayed to FCI S.A.S. and the device should not be used if the patient's safety is threatened.

#### Precautions for use

Do not use excessive force or bend the metal guide-tube when using the LACRIJET®, as this can weaken the guide-tube, increasing the danger of breakage, and makes the ejection of the stent difficult. The device should be stored between 10°C/50°F and 30°C/86°F and away from sunlight. Keep dry. Remove the device from its packaging in aseptic conditions.

#### Adverse effects

As in any type of surgery, there are risks related to the material and/or to developments of the initial pathology. Potential complications associated with the implantation of the device include, but are not limited to the following:

Intraoperative complications:

- False passage
- Epistaxis
- Excessive punctum dilation
- Stent dislocation from inserter

Postoperative complications:

- Exteriorization or loss of the tube,
- Canalicular migration of the tube,

- Keratopathy or ulceration,
- Corneal erosion,
- Canaliculitis,
- Foreign body sensation (irritation, pruritus...)
- Granuloma.

Complications can be handled with appropriate care.

The use of eye lotions containing methylene blue may cause coloration of the stent. This does not present a risk for the tolerance of the stent.

Any complication that occurs in relation with the device should be notified to FCI S.A.S. and the competent authority of the Member State in which the user and/or the patient is established, or, for users and/or patients not subject to EU requirements, to FCI S.A.S. and to the applicable regulatory authority in conformity with the requirements of that regulatory authority.

### Instructions for use

#### Exploration:

As with any other nasocanalicular intubation, the lacrimal ducts must be explored first with a surgical instrument as a Bowman probe before inserting the LACRIJET®. The result of this exploration will determine which stent and of what length of stent is suitable for the patient. The location and the severity of the obstructions are identified.

During the exploration, it is important to avoid false passages and submucosal passages.

Proper placement will be checked by carefully introducing a wide metallic probe with a blunt tip under lower nasal meatus in search of the metal-to-metal contact of the 2 probes (Fig. 1).

#### Selecting the model:

In order to determine which LACRIJET® model (30-50mm) is suitable for the treatment of stenosis, the distance between the lacrimal punctum and the nasal floor is measured. The length of the stent should be equal or superior to this distance (Fig.3).

SIZER FOR LACRIJET® (ref S1.1288) is used to select the appropriate length:

Insert the measuring device until it reaches the nasal cavity floor. Confirm the appropriate placement with another blunt-tipped metallic probe by achieving metal-to-metal contact with the measuring device (Fig.1).

Locate the marking on the measuring device that corresponds to the punctal opening.

Select the appropriate LACRIJET® model according to the markings on the measuring device (Fig.3).

Note: LACRIJET® 15 and 20mm are primarily used for trauma cases. Selection is based on the length of the laceration and lacrimal anatomy. The choice of 15 or 20mm length is decided on a case-by-case basis, using the surgeon's best clinical judgement.

#### Insertion:

##### **Lacrimal duct obstruction (30 mm – 50 mm sizes)**

A/ Lacrimal catheterization with the fitting guide

- Dilation of the lacrimal punctum,

- Search for the contact with the bone, then rotate the guide and proceed the nasolacrimal catheterization until reaching the nasal cavity floor. Proper placement will be checked by carefully introducing a wide metallic probe with a blunt tip under lower nasal meatus in search of the metal-to-metal contact of the 2 probes (Fig.1). The guide is then pulled back for one millimeter

B/ Ejection of the LACRIJET®

During the whole procedure, the body of handpiece should be hold firmly and stay in the axe of lacrimal sac.

Based on the selection of the LACRIJET® length, one of the two scenarios should be followed:

- Scenario 1: When the length from the lacrimal punctum to the nasal floor is equal to the length of the LACRIJET® and the fixation head is in contact with the lacrimal punctum (Fig.8):
  1. Removal of the unlocking part (F)
  2. Release the piston.
  3. Keep pressure on the fitting guide inside the lacrimal duct, to so the fixation head remains in contact with the lacrimal punctum
  4. Continue the piston's complete release.
- Scenario 2: When the length of the LACRIJET is greater than the distance between lacrimal punctum to nasal floor, and the fixation head is not in contact with the lacrimal punctum (Fig.9):
  1. Remove the unlocking part (F)
  2. Pull the piston (G) approximately 5 mm. This will retract the guide into the body of the handpiece and start the ejection of the LACRIJET®' stent. It is important to begin retraction of the guide to avoid pushing the device into the nasal floor.
  3. While releasing the piston, slowly push the fitting guide inside the lacrimal duct, and stop when the fixation head is in contact with the lacrimal punctum.
  4. Continue to completely release the piston and remove the handpiece.
  5. Apply the fixation head against the lacrimal punctum using the tip of FCI plug inserter (ref. S1.3090), then completely remove the handpiece.

Important: the fixation head must be in contact with the lacrimal punctum, without tending to spontaneous migration (risk of secondary intra-lacrimal migration), or to remain in distance (risk of premature expulsion if the LACRIJET® is pressed while in the canaliculus). In both cases, it is recommended to withdraw the LACRIJET® and restart the procedure from the beginning.

C/ Punctal fixation of the LACRIJET®

- The metallic tip of the plug inserter is introduced in the lumen if the fixation head.

- A slight back-pull on the eyelid while the bulb is pushed (as a Monoka® or Masterka® type stent) inside the vertical canaliculus.

- Pressure applied on the two wings of the plug inserter will disengage the fixation head from the plus inserter.

Important: ensure that the collarette is evenly applied on the punctal ring. The shape of the fixation head keeps the probe from migrating and from expulsing.

##### **Repair of Canalicular Lacerations (15 and 20 mm sizes):**

The punctum is dilated, with a tapered probe, to a sufficient diameter to facilitate insertion of the LACRIJET®. The stent should be ready for insertion since the punctum will rapidly constrict after dilation.

LACRIJET® is used like a metal probe to facilitate the catheterization between the distal and proximal portions of the canaliculus. Once the probe has been fully advanced through the canaliculus, the position of the plug should be immediately above the punctal opening. To release the stent:

1. Remove the unlocking part (F)
2. Release the piston.
3. Continue to completely release the piston and remove the handpiece.
4. Apply the fixation head against the lacrimal punctum using the tip of FCI plug inserter (ref. S1.3090), then completely remove the handpiece.

Remove of the LACRIJET®:

- No anesthesia is necessary.

- The collarette that emerges from the punctum is grasped with forceps.

- A traction removes the LACRIJET® from the lacrimal duct.

- An inspection ensures that the stent is completely removed.

#### ATTENTION: TECHNICAL MISTAKES:

1. False Passage: (Fig.1, Fig.5)

The detection of false passage during the instrument exploration step and during the use of the measuring device or after insertion of the LACRIJET® is realized by searching metal-to-metal contact between the probe, the sizer or LACRIJET® guide tube and a second blunt-tipped probe inserted under the lower nasal meatus.

2. In the case where the fixation head is not in contact with the punctum after step 4 (Fig. 6)

It is because the fitting guide, which was just pushed, comes up against the nasal cavity floor. Put once again the piston for approximately 5mm, then push the fitting guide inside the lacrimal duct. Then, with one or two attempts, the fixation head should be able to get in contact with the punctum.

3. Wrong Position: (Fig.7)

The stent is totally released from the handpiece and the fixation head remains 2 to 3 mm above the punctum. One should not try to force the stent in further with the applicator, instead verify the length of the LACRIJET® stent and restart the procedure again.

### Performances, mode of action and clinical benefits

LACRIJET is used for monocular intubation. The silicone stent acts as a conformer. In cases of canalicular lacerations the silicone tube guides wound healing and reduces post-traumatic shrinkage. The shape of the fixation head keeps the probe from migrating and expulsing. The Summary of Safety and Clinical Performance is available via Eudamed.

### Compatibility with other devices

The implanted part of LACRIJET® (stent) is MRI compatible.

### Information to be communicated to the patient

Do not rub the eye with LACRIJET® stent. An implant card is provided in the box of the product. It is intended to be filled out by the healthcare providers and given to the patient. It allows the traceability of the implanted product. Any serious incident that occurs in relation with the device should be notified to FCI S.A.S. and the competent authority of the Member State in which the patient is established. The user must inform the patient of the contraindications and potential complications related to the device that are set out in the present document.

### Waste disposal

With household waste or recycle	With hospital waste (contaminated products)
Outer box, instructions for use and traceability labels	Stent
Blisters	Handpiece

The disposal of this device does not present any physical risk.

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©LACRIJET