fr	NOTICE D'UTILISATION – KIT BILATÉRAL OPHTACATH®	2
en	INSTRUCTION FOR USE – BILATERAL OPHTACATH® KIT	3
de	GEBRAUCHSANWEISUNG – BILATERALES OPHTACATH® SET	5
it	ISTRUZIONI PER L'USO – KIT BILATERALE OPHTACATH®	7
es	INSTRUCCIONES DE USO – KIT BILATERAL OPHTACATH®	9
pt	INSTRUÇÕES DE USO – KIT BILATERAL OPHTACATH®	11
nl	GEBRUIKSAANWIJZING – BILATERALE OPHTACATH® KIT	13
da	BRUGSVEJLEDNING – BILATERALT OPHTACATH® SÆT	15
fi	KÄYTTÖOHJEET – BILATERAL OPHTACATH® -SETTI	16
sv	BRUKSANVISNING – BILATERAL OPHTACATH® SATS	18
el	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ – ΑΜΦΙΠΛΕΥΡΟ ΚΙΤ ΟΡΗΤΑCATH®	20
et	KASUTUSJUHEND – BILATERAALNE KOMPLEKT OPHTACATH®	22
lv	LIETOSANAS INSTRUKCIJA – BILATERAL OPHTACATH® KOMPLEKTS	24
lt	NAUDOJIMO INSTRUKCIJOS – DVIPUSIS KOMPLEKTAS „OPHTACATH®“	25
pl	INSTRUKCJA UŻYTKOWANIA – DWUSTRONNY ZESTAW OPHTACATH®	27
sl	NAVODILA ZA UPORABO – DVOSTRANSKI KOMPLET OPHTACATH®	29
sk	NAVOD NA POUZIVANIE – OBOJSTRANNÁ SÚPRAVA OPHTACATH®	31
cs	NAVOD K POUZITÍ – OBOUSTRANNÁ SADA OPHTACATH®	33
ro	INSTRUCȚIUNI DE UTILIZARE – KIT BILATERAL OPHTACATH®	34
bg	ИНСТРУКЦИИ ЗА УПОТРЕБА – ДВУСТРАНЕН КОМПЛЕКТ ОΡΗΤΑCATH®	36
hr	UPUTE ZA UPOTREBU – KOMPLET BILATERAL OPHTACATH®	38
hu	HASZNALATI UTASITAS – KÉTOLDALI OPHTACATH® SZETT	40
no	BRUKSANVISNING – BILATERALT OPHTACATH®-SETT	42
tr	KULLANIM – BİLATERAL OPHTACATH® KİTİ	43
zh	使用说明 - OPHTACATH® 双侧套件	45
ar	تعليمات الاستخدام – مجموعة OPHTACATH® الثنائية	
Illustration.....		45
Symbols.....		49
IFU-2138 313 46_C		

FCI S.A.S.
20-22, rue Louis Armand
75015 PARIS – FRANCE

☎ : + 33 (0) 1 53 98 98 98

📠 : + 33 (0) 1 53 98 98 99

@ : fci@fci.fr

🌐 : www.fci.fr

Distributed in U.S.A. by:

FCI Ophthalmics Inc.
30 Corporate Park Drive
Suite # 310/320
Pembroke, MA 02359

@ : info@fci-ophthalmics.com

🌐 : www.fci-ophthalmics.com

☎ : 800-932-4202

📠 : 781-826-9062

R_x Only

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en

Instructions for use: 2138 313 46_C

Date of revision of the instructions for use : 2023-11

BILATERAL OPHTACATH® KIT

Description and packaging

Product scope of this IFU:

- S1.4112: OPHTACATH® 2mm STERILE BALLON CATHETER
- S1.4113: OPHTACATH® 3mm STERILE BALLON CATHETER
- S1.4122 : OPHTACATH® 2mm Bilateral Kit.
- S1.4132: OPHTACATH® 3mm Bilateral Kit.

OPHTACATH® is a sterile and single-use balloon catheter consisting of a semi-flexible stainless steel stylet covered by a PBX tubing that ends with a nylon balloon. The balloon is covered by a protective polyethylene sleeve. The balloon is designed to be inflated to a known diameter and length at the specified pressure. Markings are placed at 10 mm and 15 mm proximal to the working portion of the balloon. The markings are placed to indicate the insertion limits of the balloon catheter. The balloon is designed to be inflated to a given diameter and length under a reference pressure of 8 atm. The overall length of the catheter is approximately 240 mm long.

OPHTACATH® 2MM STERILE BALLON CATHETER		
	Before inflation	During inflation
Diameter	0.90 mm	2 mm
Length	13 mm	13 mm

OPHTACATH® 3MM STERILE BALLON CATHETER		
	Before inflation	During inflation

Diameter	1.0 mm	1.0 mm
Length	15 mm	15 mm

The device is sterilized using ethylene oxide in its final simple pouch packaging. A green indicator shows that the product followed a sterilization cycle validated by FCI. The device is supplied by box of 1. The device is for single use only and is not resterilizable.

Intended purpose

The OPHTACATH® kits are intended for use during dilation of the obstructed nasolacrimal ducts to treat epiphora.

Indications

Medical indication :

Nasolacrimal duct obstruction

Intended patient population:

Children (12 months and older) and adults.

Body parts in contact with the device:

Nasolacrimal ducts.

Intended users:

Surgeon with special training for nasolacrimal duct surgery.

Intended use environment :

Operating room

Contraindications

Utilization of the device is contraindicated in the following cases:

- Lacrimal sac tumor
- Dacryoliths

Warnings

- Do not touch the balloon with instruments which may puncture the balloon.
- The device is for single use and must not be resterilized. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety.
- The balloon must not be removed while inflated.
- Do not inflate the balloon above 15 atm (15 bars). Inflation of the balloon catheter above 15 atm (15 bars) may cause rupture of the balloon tip and result in serious injury.
- Do not forcibly push on the control button to lock or unlock the inflation device. Doing so may irreversibly damage the plunger, and the inflation device may not perform as intended.
- Check packaging integrity before use, and do not use the device in the event of packaging being damaged, unintentionally opened before use, or if it was exposed to environmental conditions outside of those specified.
- Do not use after the expiration date shown on the package.
- The guaranteed maximum residual rate of ethylene oxide at the time the device was placed on the market was 60 µg/MD per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

Precautions for use

The device should be stored between 5°C/41°F and 45°C/113°F and away from sunlight. Keep dry.

Remove the device from its packaging in aseptic conditions.

Adverse effects

As in any type of surgery, there are risks linked to the material or to developments of the initial pathology.

Potential complications associated with the use of the device include, but are not limited to the following:

- Swelling.
- Epistaxis
- Acute dacryocystitis

Any complication that occurs in relation with the device should be notified to FCI S.A.S. and the competent authority of the Member State in which the user and/or the patient is established, or, for users and/or patients not subject to EU requirements, to FCI S.A.S. and to the applicable regulatory authority in conformity with the requirements of that regulatory authority.

Instructions for use

1. Select size of balloon catheter that is appropriate for the patient. Before insertion of the balloon catheter, remove the protective sleeve from the distal end of the balloon catheter.
2. Lubricate the balloon catheter with ophthalmic ointment if needed.
3. Before each inflation, the manometer of the inflation device must be set at 0 ± 1 atm (bars) of pressure
4. Under anesthesia, dilate the punctum, as necessary. Pass an appropriately sized Bowman probe into the nasolacrimal duct. Confirm the presence of the probe in the nose. Remove the probe.
5. Pass the balloon catheter through the canaliculus and lacrimal sac into nasolacrimal duct until the nasal floor is reached. Confirm the presence of the distal part of the balloon catheter in the nose. Pull the balloon proximally and align the punctum with the marking located at 15mm proximal to the balloon (Figure A).
6. Prepare sterile water or saline following the manufacturer's recommendations.
7. Press button to unlock the inflation device; (Figure D).
8. Hold the inflation device with its manometer downward and fill it with 10 cc of sterile water or saline.
9. Hold the inflation device with its manometer upward to purge the air contained inside.
10. Release button to unlock the screw rod (Figure D).
11. Connect the filled inflation device to the 4 way stopcock and the balloon catheters (Figure C).
12. Inflate the balloon gradually by slowly turning the handle clockwise up to a fill pressure of 8 atm (bars) and hold for 90 seconds (Figure C).
13. Deflate the balloon by turning the handle counterclockwise until the inflator gauge of the manometer is in the red area (0 atm).
14. If necessary, an additional inflation can be done by gradually inflating the balloon to 8 atm (bars) and holding for 60 seconds, and then deflating the balloon as in Step 13.
15. Pull the balloon proximally and align the punctum with the marking located at 10mm proximal to the balloon catheter (Figure B). The balloon is gradually inflated to 8 atm (bars) and held for 60 seconds.
16. Deflate the balloon in the same way as Step 13.
17. If necessary, an additional inflation can be done by gradually re-inflating the balloon to 8 atm (bars) and holding for 40 seconds.

18. After the last inflation is performed, fully deflate the balloon by turning the handle counterclockwise until the inflator gauge of the manometer is in the red area (0 atm) and all fluid has been aspirated out of the balloon.

19. Withdraw the catheter progressively.

Performances, mode of action and clinical benefits

The clinical benefit provided by OPHTACATH® is the improvement of the epiphora via the restoration of the flow of tears. The mode of action of the device is mechanical: the catheter is inflated and deflated at the level of the obstruction in order to dilate the lacrimal ducts.

Compatibility with other devices

OPHTACATH® Balloon Catheter is intended to be used with the accessories provided with it in the kit, as following:

- The bilateral 4-way stopcock (Merit Medical) is intended to be used as an accessory within fluid management.
- Balloon Inflation Device is intended for use during nasolacrimal duct dilatation procedures to inflate balloon catheters.

Information to be communicated to the patient

Any serious incident that occurs in relation with the device should be notified to FCI S.A.S. and the competent authority of the State in which the patient is established. The user must inform the patient of the contraindications and potential complications related to the device that are set out in the present document.

Waste disposal

Used devices may be a potential biohazard. Handle and discard used devices in accordance with the accepted medical practice and applicable local, state, and federal laws and regulations.

With household waste or recycle	With hospital waste (contaminated products)
Outer box, instructions for use and traceability labels	OPHTACATH® balloon catheter
Pouches	The bilateral 4-way stopcock
	Balloon inflation device

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