REF. 2138 301 02 M

FR	NOTICE D'UTILISATION - PERFLUOROCARBONES LIQUIDES	2
EN	INSTRUCTION FOR USE - PERFLUOROCARBON FLUIDS	4
DE	GEBRAUCHSANWEISUNG - FLÜSSIGE PERFLUOROKARBONE	5
IT	ISTRUZIONI PER L'USO - PERFLUOROCARBONI LIQUIDI	7
ES	INSTRUCCIONES DE USO - PERFLUOROCARBONOS LÍQUIDOS	9
PT	INSTRUÇÕES DE USO - PERFLUOROCARBONOS LÍQUIDOS	11
NL	GEBRUIKSAANWIJZING - VLOEIBARE PERFLUOROCARBONATEN (PFCs)	13
DA	BRUGSVEJLEDNING - PERFLUOROCARBON-VÆSKER	15
FI	KÄYTTÖOHJEET - NESTEMÄISET PERFLUORATUT HIILIVEDYT	17
SV	BRUKSANVISNING - FLYTANDE PERFLOURKARBON	19
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - ΥΓΡΑ ΥΠΕΡΦΘΟΡΑΝΘΡΑΚΑ	21
LV	LIETOŠANAS INSTRUKCIJA - PERFLUOROGLEKĻA ŠĶIDRUMI	23
PL	INSTRUKCJA UŻYTKOWANIA - PŁYNY PERFLUOROWĘGLOWE	25
SL	NAVODILA ZA UPORABO - PERFLUOROOLJIKOVODIKOVE TEKOČINE	27
SK	NÁVOD NA POUŽITIE – PERFLUÓROKARBÓNÓVÉ KVAPALINY	29
CS	NÁVOD K POUŽITÍ – PERFLUOROUHLÍKOVÉ KAPALINY	31
HU	HASZNÁLATI ÚTMUTATÓ - PERFLUORKARBON FOLYADÉKOK	33
TR	KULLANIM - PERFLUOROKARBON SIVILAR	35
AR	إرشادات الاستخدام - سوائل بيرفلوروكربون	36
ZH-CN	使用说明 - 眼科手术用重水 (FCI-DECA) 使用说明 - 眼科手术用重水 (FCI-OCTA)	38
ZH-TW	使用說明 - 全氟碳化物液體 (PFCL)	40

Symbols..... 43

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PERFLUOROCARBON FLUIDS FCI-DECA and FCI-OCTA are adjuvants in vitreo-retinal surgery. They are available in two versions:

FCI-DECA: purified perfluorodecalin, available in different capacity vials (5ml, 7ml).

FCI-OCTA: purified perfluoro-N-octane, available in different capacity vials (5ml, 7ml).

PERFLUOROCARBON FLUIDS FCI-DECA and FCI-OCTA are supplied purified, sterile and non-pyrogenic. This product does not contain latex nor phthalates. They are packaged in glass vials, sealed with a polymer stopper and a metal crimp cap. They are sterilized in their final package to guarantee sterility of the contents.

The PFCL FCI-DECA and FCI-OCTA are subjected to a sterilizing filtration and steam sterilization in order to facilitate handling under aseptic conditions.

PHYSICAL PROPERTIES

	Refractive index	Boiling point	Vapor tension	Specific gravity
FCI-DECA	1.31	141 °C	13.5 mmHg / 37 °C	1.94 g/cm ³
FCI-OCTA	1.27	102 °C	57 mmHg / 37 °C	1.76 g/cm ³

Intended purpose:

PFCL (FCI-DECA and FCI-OCTA) are used as per-operative devices for transitory intra-ocular endotamponade.

Indications :

Medical conditions: PERFLUOROCARBON FLUIDS FCI-DECA and FCI-OCTA are used as per-operative devices, for transitory intra-ocular endotamponade; in particular with giant tear, retinal detachment (RD), with proliferative vitreal retinopathy (PVR) including diabetic proliferative retinopathy, flushing out of the subretinal liquids, retrieval of foreign bodies, and retrieval of a luxated natural or artificial lens in the vitreous.

Intended patient population: children (12 months and older) and adults.

Body parts in contact with the device: the retina and sclera.

Intended users: surgeons trained in retinal surgery.

Intended use environment: operating room

Contraindications :

Utilization of the device is contraindicated in the following cases:

The product should not be allowed to infiltrate under sub-retinal tissues since there are no specific studies on tolerance. Retina too fibrotic to enable reattachment. Posterior dehiscence.

Warnings:

The PFCL should not be injected directly into the vortex vein. In addition, inadvertent leakage of intraocular perfluorocarbon fluids into the vortex vein should be avoided.

The device is for single use and must not be resterilized.

Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety

Precautions for use:

The PFCL FCI-DECA and FCI-OCTA are restricted for use by surgeons trained in retinal surgery.

The PFCL FCI-DECA and FCI-OCTA are per-operative devices for vitreo-retinal surgery.

The volume injected must be adapted to each procedure.

PFCL FCI-DECA and FCI-OCTA are designed to be removed at the end of the procedure.

Do not use a glass syringe under any conditions.

Do not use the needle (21G) provided with PFCL for intraocular injection of the product.

The PFCL significantly decrease ultrasound speed.

Before use, check the individual package to be sure no seals have been broken and the sterility is still valid.

PERFLUOROCARBON FLUIDS FCI-DECA and FCI-OCTA are single-use liquids and must not be re-sterilized. The sterility and purity of the product are only guaranteed for a single use. If you have some product left, do not use it for another patient. They must not be used after the expiration date shown on the package. No particular storage requirement. Store at room temperature.

Special precautions are indicated when PFCL may enter the vascular circulation (e.g., during endoresection of intraocular tumors and ocular trauma).

Adverse effects :

As in any type of surgery, there are risks linked to the material or to developments of the initial pathology. Potential complications associated with the use of the device include, but are not limited to the following:

- Cataract
- IOP raised
- Low pressure ocular
- PFCL in the anterior chamber
- Anterior uveitis
- Posterior segment PFCL
- Subretinal PFCL
- Cystoid macular oedema
- ERM
- Redetachment
- Intraocular lens repositioning
- Glaucoma
- Bullous kerathopathy
- Hyphema
- Vitreous hemorrhage
- Blurred/Loss vision
- Full macular hole

- Inflammation
- Eye injury
- Discomfort
- Toxic Anterior Segment Syndrome (TASS)
- Air embolism

Any complication that occurs in relation with the device should be notified to FCI S.A.S. and the competent authority of the Member State in which the user and/or the patient is established, or, for users and/or patients not subject to EU requirements, to FCI S.A.S. and to the applicable regulatory authority in conformity with the requirements of that regulatory authority.

Instructions for use :

Fit the needle provided onto the syringe and pierce the stopper after having removed the plastic part of the metal cap of the vial and pierce the stopper. Transfer the product from the vial into the syringe. Do not use this needle for intra-ocular injection of the product: Replace the needle with a blunt cannula for intra-vitreous injection.

- Giant tear, retinal detachment (RD)

Inject PFCL FCI-DECA and FCI-OCTA slowly and progressively, just in front of the tear, after total pars plana vitrectomy. Maintain the tip of the cannula inside of the bubble during the injection to prevent the product dispersion. If despite this precaution, dispersion occurs, the droplets will spontaneously regroup into the principal PFCL FCI-DECA and FCI-OCTA bubble or can be removed thanks to an extrusion cannula. Thanks to their high density PFCL FCI-DECA and FCI-OCTA reattach the retinal flap from the posterior pole towards the periphery, flushing ahead the sub-retinal liquid, the removal of which is then passively performed through the sclerotomy incisions.

The progressive retina re-attachment should be controlled to avoid sub-retinal passage: any residual traction should be cancelled in order to get a soft re-attachable retina. Always keep PFCL FCI-DECA and FCI-OCTA in front of the edge of the tear to prevent sub-retinal passage.

- RD with vitreo-retinal proliferation (PVR)

After vitrectomy, proceed with the dissection of the posterior epiretinal membranes. Then, inject PFCL FCI-DECA and FCI-OCTA slowly and gradually, restrain to completely fill the vitreous cavity. The posterior retina is stabilized and allows the treatment of the anterior epiretinal membranes. Always proceed ahead of the bubble to avoid the sub-retinal passage in case of iatrogenic tear occurrence.

In the case of tight retinal funnel, with strong anterior pull, open the funnel with an injection of a small volume of PFCL FCI-DECA and FCI-OCTA : the surface of the retina is stabilized and enables the anterior epiretinal membranes dissection.

Then, after removal of the product, proceed with the posterior membranes dissection. When the retina is completely released, a second injection to totally fill the cavity enables the re-attachment and the sub-retinal fluids drainage as in the giant tear case.

- Retinopexy using endophotocoagulation can be performed

The contact of the product with the retina should be limited as much as possible. It should not be left in place in the eye at the end of the procedure. The total quantity injected shall be removed at the end of the surgery.

The removal is carried out, either passively through the downwards fluid, air, perfluorinated gas or silicone oil injection: PFCL FCI-DECA and FCI-OCTA is thus pushed downwards and passively removed through an extrusion cannula, or by aspiration.

A meniscus may occur at the end of the PFCL FCI-DECA and FCI-OCTA -air or perfluorinated gas exchange: the reinjection of a small volume of Ringer or BSS disperses the residual quantity and enables the complete removal.

PFCL FCI-DECA and FCI-OCTA -air or perfluorinated gas exchange duration should not be extended; high intraocular pressure should be maintained during the exchange.

The risk of sub-retinal passage is increased in the case of wide peripheral retinotomy, especially when it reaches 360 degrees. Rinsing with a saline solution during the exchange seems to be useful to reduce the sub-retinal passage.

Foreign bodies and subluxated natural crystalline or artificial lens can be removed thanks to PFCL FCI-DECA and FCI-OCTA once they are completely released from surrounding vitreous after performing a complete pars plana vitrectomy. The PFCL FCI-DECA and FCI-OCTA is directly injected underneath the material to be removed, which then will float on PFCL FCI-DECA and FCI-OCTA surface and can then be safely extracted.

Performance, mode of action and clinical benefits

PFCL is designed to reapply the retina by mechanical pressure against the choroid based on its high density.

Compatibility with other devices

Compatibility with MRI imaging: The PFCL (FCI-DECA and FCI-OCTA) are perfluorinated compounds and are characterized by the lack of hydrogen. The PFCL do not contain any hydrogen atoms and therefore do not give any signal.

The products are provided with a single needle (21G) and syringe (5 ml or 10 ml) single use.

Information to be communicated to the patient:

The patient should be informed of possible residual droplets and their eventual side effects that they should bring to the caretaker's attention.

Any serious incident that occurs in relation with the device should be notified to FCI S.A.S. and the competent authority of the State in which the patient is established.

The user must inform the patient of the contraindications and potential complications related to the device that are set out in the present document.

Waste disposal

With household waste or recycle	With hospital waste (contaminated products)
Outer box, instructions for use and traceability labels	Retreated solution
Pouches	Vial with remaining solution
	Needle
	Syringe with remaining solution

The disposal of this device does not present any physical risk.

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