FR	NOTICE D'UTILISATION – PROTECTEUR OCULAIRE	2
EN	INSTRUCTION FOR USE – OCULAR SHIELD	3
DE	GEBRAUCHSANWEISUNG – OKULARSCHUTZKAPPE	4
IT	ISTRUZIONI PER L'USO – SCUDO OCULARE	5
ES	INSTRUCCIONES DE USO – PROTECTOR OCULAR	6
PT	INSTRUÇÕES DE USO – PROTETOR OCULAR	7
NL	GEBRUIKSAANWIJZING – OOGBESCHERMER	8
DA	BRUGSVEJLEDNING – ØJENSKJOLD	9
FI	KÄYTTÖOHJEET – SILMÄSUOJUS	10
SV	BRUKSANVISNING – OKULÄR SKÖLD	11
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ – ΕΝΔΟΦΘΑΛΜΙΟ ΚΑΛΥΜΜΑ	12
ET	KASUTUSJUHEND – SILMAKILP	13
LV	LIETOSANAS INSTRUKCIJA – ACS AIZSARGS	14
LT	NAUDOJIMO INSTRUKCIJOS – AKIES APSAUGOS PRIEMONĖ	15
PL	INSTRUKCJA UŻYTKOWANIA – OSŁONA OCZU	16
SL	NAVODILA ZA UPORABO – ŠČITNIK ZA OČI	17
SK	NAVOD NA POUZIVANIE – OČNÝ ŠTÍT	18
CS	NAVOD K POUZITI – OČNÍ ŠTÍT	19
RO	INSTRUCȚIUNI DE UTILIZARE – PROTECTOR OCULAR	20
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА – ОЧЕН ПРОТЕКТОР	21
HR	UPUTE ZA UPOTREBU – OČNI ŠTIT	22
HU	HASZNALATI UTASITAS – SZEMVÉDŐ	23
NO	BRUKSANVISNING – ØYESKJOLD	24
TR	KULLANIM – OCULAR SHIELD	25
ZH	使用说明书 – 护眼罩 (OCULAR SHIELD)	26
Illustrations.....		27
Symbols.....		28
<b>REF. 2138 311 18 C</b>		

FCI S.A.S.  
20-22, rue Louis Armand  
75015 PARIS – FRANCE

☎ : + 33 (0) 1 53 98 98 98

📠 : + 33 (0) 1 53 98 98 99

@ : fci@fci.fr

🌐 : www.fci.fr

Distributed in U.S.A. by:

FCI Ophthalmics Inc.  
30 Corporate Park Drive  
Suite # 310/320  
Pembroke, MA 02359

☎ : 800-932-4202

📠 : 781-826-9062

@ : info@fci-ophthalmics.com

🌐 : www.fci-ophthalmics.com

Rx Only

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The device is for single use only.

- S6.2301: OCULAR SHIELD

OCULAR SHIELD is a protective egg-shaped shield (maximal diameter: 27mm) made of PMMA (Polymethylmetacrylate), with a prehension pin which is used to place and remove the shield easily.

The device is sterilized using Gamma radiation in its final simple blister packaging. A red indicator shows that the product followed a sterilization cycle validated by FCI. The device is supplied by box of 1. The device is for single use only and is not resterilizable.

### Indications

OCULAR SHIELD is designed to protect the eyeball during conventional eyelid surgery or lacrimal duct surgery (DCR).

This medical device is for professional use only.

### Contraindications

No known contraindications.

### Warnings

The device is for single use and must not be resterilized. Reuse and/or reesterilization may compromise device performance, which could cause serious harm to the patient's health and safety.

Check packaging integrity before use, and do not use the device if package is damaged.

Do not use the product if the indicator is not red, as it could endanger the patient.

Do not use after the expiration date shown on the package.

The device is intended for use for not more than 30 days.

In case the user identifies faulty labelling information, it shall be relayed to FCI S.A.S. and the device should not be used if the patient's safety is threatened.

### Precautions for use

The device should be stored between 10 and 30°C and away from sunlight. Keep dry. Remove the device from its packaging in aseptic conditions.

### Adverse effects

As in any type of surgery, there are risks linked to the material or to developments of the initial pathology. Potential complications associated with the use of the device include, but are not limited to the following:

- Corneal abrasion

- Corneal infection

- Discomfort

- Red eyes and photophobia

- Sub conjunctival hemorrhage

Any complication that occurs in relation with the device should be notified to FCI S.A.S. and the competent authority of the Member State in which the user and/or the patient is established.

### Instructions for use

Use the OCULAR SHIELD prehension pin to place and remove it easily.

### Performances, mode of action and clinical benefits

OCULAR SHIELD forms a mechanical barrier between the eyelids and the cornea, which allows the surgeon to protect the eyeball of the patient during the surgery.

### Compatibility with other devices

FCI S.A.S. does not claim the compatibility of OCULAR SHIELD with other devices.

### Information to be communicated to the patient

Any serious incident that occurs in relation with the device should be notified to FCI S.A.S. and the competent authority of the Member State in which the patient is established. The user must inform the patient of the contraindications and potential complications related to the device that are set out in the present document.

### Waste disposal

With household waste or recycle	With hospital waste (contaminated products)
Outer box, instructions for use	Ocular shield
Blister	

The disposal of this device does not present any physical risk.

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