fr	NOTICE D'UTILISATION – FCI Protect®	2
en	INSTRUCTION FOR USE – FCI Protect®	3
de	GEBRAUCHSANWEISUNG – FCI Protect®	4
it	ISTRUZIONI PER L'USO – FCI Protect®	6
es	INSTRUCCIONES DE USO – FCI Protect®	7
pt	INSTRUÇÕES DE USO – FCI Protect®	8
nl	GEBRUIKSAANWIJZING – FCI Protect®	10
da	BRUGSVEJLEDNING – FCI Protect®	11
fi	KÄYTTÖOHJEET – FCI Protect®	12
sv	BRUKSANVISNING – FCI Protect®	14
el	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ – FCI Protect®	15
et	KASUTUSJUHEND – FCI Protect®	16
lv	LIETOSANAS INSTRUKCIJA – FCI Protect®	18
lt	NAUDOJIMO INSTRUKCIJOS – FCI Protect®	19
pl	INSTRUKCJA UŻYTKOWANIA – FCI Protect®	20
sl	NAVODILA ZA UPORABO – FCI Protect®	21
sk	NAVOD NA POUZIVANIE – FCI Protect®	23
cs	NAVOD K POUZITÍ – FCI Protect®	24
ro	INSTRUCȚIUNI DE UTILIZARE – FCI Protect®	25
bg	ИНСТРУКЦИИ ЗА УПОТРЕБА – FCI Protect®	27
hr	UPUTE ZA UPOTREBU – FCI Protect®	28
hu	HASZNALATI UTASITAS – FCI Protect®	29
no	BRUKSANVISNING – FCI Protect®	31
tr	KULLANIM – FCI Protect®	32
zh	使用说明- FCI 保护® FCI Protect®	33
Illustrations.....		35
Symbols.....		39
REF. 2138 315 15 C		

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en

Instructions for use: 2138 315 15 C

Date of revision of the instructions for use: 2022-05

FCI PROTECT®

Description and packaging:

Product scope of current IFU:

Product Code	Designation
S5.9100	FCI Protect® 2%

FCI Protect® is a sterile syringe for single use containing 2 ml of solution. It's a solution isotonic, clear, viscoelastic and sterile.

1 ml of solution contains 20mg HPMC (hydroxypropyl methylcellulose) as rheological active ingredient dissolved in balanced salt solution containing sodium chloride, potassium chloride, calcium chloride, magnesium chloride, sodium citrate, sodium acetate, and water for injection. The package includes a cannula 23G.

Solution's Characteristics

Rheologically active component	Hydroxypropyl methylcellulose (HPMC)
Concentration	2.0% HPMC
Volume	2.0 ml
pH	6.8-7.6
Average molecular weight Mw (Dalton)	86 000
Polydispersity	1-5
Osmolality (mOsmol/kg)	200 – 400

FCI Protect® is sterilized by moist heat inside its sterile packaging of single blister.

Intended purpose:

FCI Protect® Syringe is intended for topical use. It serves as an epithelium protection by coating and lubrication of the eye during the use, for example, of vitrectomy lenses, or during gonioscopy exam (exam of the fundus of the eye and the iridicorneal angle with three-mirror contact lens).

Indications:

FCI Protect® is indicated for the protection of the cornea during the use of vitrectomy lenses or during gonioscopy exam.

Contra-indications:

FCI Protect® should not be used in cases of a hyper-sensitivity against any components in the material.

Warnings

FCI Protect® is a single use device that must not be re-sterilized. Re-use and/or re-sterilization could affect the device performances and compromise the patient's safety and health.

During the surgery, a prolonged exposure of FCI Protect® should be evited, it may cause the deleterious effects upon the corneal endothelium.

Failure to follow the assembly instructions, or use of an alternate cannula may result in cannula detachment and the possibility of serious injury of patient.

In case the user identifies faulty labelling information, it shall be relayed to FCI S.A.S. and the device should not be used if the patient's safety is threatened.

Check packaging integrity before use, and do not use the device in the event of packaging being damaged.

Do not use FCI Protect® after the expiration date as it might compromise the product's efficiency and the patient's safety.

Precautions for use

This medical device is intended for use by ophthalmologists only. FCI Protect® is provided with an optional cannula.

The device should be stored between 2 -25 °C (36-77°F) in a cool, dry place and away from sunlight.

Adverse effects

The other adverse reactions include: corneal oedema, Inflammation or infection.

Any complication that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Instructions for use :

Assembly

Figure 1: Hold the Luer-Lock adapter as shown in (1). Twist the cap carefully with the other hand in an anti-clockwise direction (2). Then remove the cap as shown in (3). Do not use a syringe if the cap is not on properly in its sterile packaging.

Hold the syringe as shown in (4). Insert the enclosed 23G cannula firmly as shown in (5) (do not use any other cannula). Hold the cannula and lock it into position by twisting slightly in anti-clock-wise direction (6).

Hold the syringe during administration as shown in Figure 2.

Injection for topical use

Inject the solution on the surface of the eye until it covers totally and evenly the cornea with the gel.

Elimination of solution for extraocular application

Flush the surface of the eye with saline solution.

Performance, mode of action and clinical benefits

FCI Protect® protects the corneal epithelium during its topical use by coating and moisturization of the eye.

Compatibility with other devices

The FCI Protect® is intended to be used with a 23G cannula for its injection.

Information to be communicated to the patient

Any serious incident that occurs in relation with the device should be reported to FCI S.A.S. and the competent authority of the Member State in which the patient is established.

Waste disposal:

Regarding the disposal of the device, see the instruction in the following table:

With household waste or recycle	With hospital waste (contaminated products)
Instructions for use and traceability label	Remaining solution + Used syringe
Blister	Removed solution

The disposal of this device does not present any particular risk.

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