Fr	NOTICE D'UTILISATION – SUBCYCLO®	2
En	INSTRUCTION FOR USE – SUBCYCLO®	4
De	GEBRAUCHSANWEISUNG – SUBCYCLO®	6
It	ISTRUZIONI PER L'USO – SUBCYCLO®	8
Es	INSTRUCCIONES DE USO – SUBCYCLO®	10
Nl	GEbruIKSAANWIJZING – SUBCYCLO®	12
Fi	KÄYTTÖOHJEET – SUBCYCLO®	14
No	BRUKSANVISNING – SUBCYCLO®	16
Tr	KULLANIM – SUBCYCLO®	18
Pt-Br	INSTRUÇÕES DE USO – SUBCYCLO®	20
Illustrations		22
Symbols		23
IFU-2138-371-40_A		

FCI S.A.S.
20-22, rue Louis Armand
75015 PARIS – FRANCE

 : + 33 (0) 1 53 98 98 98

 : + 33 (0) 1 53 98 98 99

 : fci@fci.fr
 www.fci.fr

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EN

Instructions for use: 2138 -371-40_A
Date of revision of the instructions for use: 2023-12

SUBCYCLO®**Description and packaging**

Product scope of this IFU: XLCYCLO: SUBCYCLO® – QUANTEL

SUBCYCLO® are composed of:

- a connector for connecting to the laser surgical machine,
- an optical fiber for laser wave transmission,
- a handle for manipulation by the surgeon,
- a removable transparent tip.

The device is sterilized using ethylene oxide in its final double pouch packaging to make it easier to handle in aseptic conditions. The device is supplied by box of 5. The device is for single use only and is not resterilizable.

Intended purpose

This product, when used with an infrared laser, is indicated for Transscleral Cyclophotocoagulation (TSCPC) of the ciliary processes.

Indications

Medical conditions: Excessive intraocular pressure (IOP).

Intended patient population: Children (12 months and older) and adults.

Body parts in contact with the device: External ocular surface (scleral conjunctiva).

Intended users: Ophthalmic surgeons.

Intended use environment: Operating room.

Contraindications

No known contraindications.

Warnings

The device is for single use and must not be resterilized. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety.

Check packaging integrity before use, and do not use the device in the event of packaging being damaged, unintentionally opened before use, or if it was exposed to storage or transport conditions outside of those specified. Do not use after the expiration date shown on the package.

If not handled properly during any step of setup (such as bending, improperly securing, poorly inserting, etc.), the glass optical fiber within the product can be damaged leading to harm to the patient or laser operator.

Excess treatment power may result in ocular surface burns or ciliary body hemorrhage. The maximum recommended laser power for this device is 2500mW. Contamination of the tip or reuse of the probe may result in burns. Heavy pigmentation at the perilimbal conjunctiva may also result in burns, therefore these areas should be avoided.

If there is an audible "pop" or "snap" during treatment, this may indicate intraocular tissue disruption; reduce power and leave the duration unchanged.

This device should be used with an appropriate laser safety filter.

No one should ever look directly into the laser light source or the tip of the probe.

All operating room personnel should be instructed to wear adapted safety goggles and avoid looking at surfaces that the laser light may be scattered from.

TSCPC treatment shall be adapted for an eye that has previously undergone surgery.

Diode laser Cyclophotocoagulation is generally not recommended in phthisical eyes with disorganized internal structures or where other alternative and more effective treatment modalities are available, especially in eyes with good vision or at risk of macular edema.

Uveitic patients or neovascular glaucoma cases must also be considered carefully due to their risk of severe inflammation and hypotony post-operatively.

The guaranteed maximum residual rate of ethylene oxide the device can release to the patient is 210µg. The residual quantity of ethylene oxide corresponds to the maximum quantity of residues that may leach over the lifetime of use.

In case the user identifies faulty labelling information, it shall be relayed to FCI S.A.S. and the device should not be used if the patient's safety is threatened.

Precautions for use

The device must be stored at temperatures between 5°C/41°F and 35°C/95°F and humidity between 30% and 75%, away from sunlight.

The device may be exposed to temperature variations between -15°C/5°F and 50°C/122°F and humidity variations between 15% and 90% during transport for 10 days.

Remove the device from its packaging in aseptic conditions.

Eye Moistness

Keep the eye moist throughout the entire laser procedure, using a viscous liquid / gel interface as needed.

Laser Settings

Recommended settings are effective for most eyes, but it is the responsibility of the clinician to select and adjust them depending on each patient characteristics and clinical history.

Cyclophotocoagulation probe inspection during procedure:

Inspect the tip of the optic fiber after each application to verify that it is not visibly blackened, is free of adherent tissue, retains its spherical shape, and is not cracked.

If the tip becomes dirty during the procedure, clean it gently with a medical alcohol swab. Solvents of flammable solutions used for cleaning should be allowed to evaporate before laser equipment is used. If discoloration on the tip cannot be removed by gentle cleaning, discard the probe. If a scleral burn occurs, discard the probe, and replace it with a new probe. Check the probe before use.

The aiming beam should emerge from the optic fiber tip with a diverging beam profile and produce a round spot on a perpendicular target with uniform brightness.

After procedure, monitoring the IOP may help prevent secondary post-operative complications.

Adverse effects

No specific adverse effects are associated with the use of SUBCYCLO®.

Several complications and side effects might be associated specifically to the surgical procedure, or to the use of LASER sources. Possible complications include but are not limited to the following:

- Conjunctival or scleral burns
- Hypotony
- Phthisis bulbi
- Scleritis
- Perforation
- Sympathetic ophthalmia
- Corneal ulcer
- Malignant glaucoma
- Scleral thinning

Any complication that occurs in relation with the device should be notified to FCI S.A.S. and the competent authority of the Member State in which the user and/or the

patient is established, or, for users and/or patients not subject to EU requirements, to FCI S.A.S. and to the applicable regulatory authority in conformity with the requirements of that regulatory authority.

Instructions for use

- Remove the probe from pouches and check it before use.
- Connect the probe to a compatible infrared laser machine (recommended wavelength: 810nm). Follow the recommendations of the laser source manufacturers (e.g.: calibration of the laser source as well as the need for protective eye wear for the corresponding wavelength of the laser model used).
- The probe is now ready to use.

1. Instructions for use of probe in continuous mode, WITH footplate attached

Laser Settings

Standard technique: Recommended settings are 1250mW for 4000 milliseconds by treatment site for a patient with a dark brown iris; and 1500mW for 2000 milliseconds by treatment site for all other iris colors. Adjust settings for every eye by increasing power as necessary.

Slow coagulation technique: Recommended settings are 1250mW for 4000 to 4500 milliseconds by treatment site for a patient with a dark brown iris; and 1500mW for 3500 to 4000 milliseconds by treatment site for all other iris colors.

Placement (Fig.1)

Hold the probe parallel to the visual axis of the eye and perpendicular to the iris plane. The shorter edge of the footplate of the probe should be held firmly against the area located between the middle and the anterior border of the limbus. Laser delivery is transscleral.

Applications and treatment (Fig.2)

Use one-half of the footplate width to determine the correct spacing between treatment sites. Align the side of the probe over the indented center of the previous treatment.

This should result in 6-7 laser applications per quadrant, 18-21 applications total, avoiding application in the 3 and 9 o'clock meridians.

2. Instructions for use of probe in pulsed mode WITHOUT footplate:

Location and marking of the ciliary body position (Fig.3)

Assess the ciliary body position thanks to the transscleral illumination technique and mark the treatment area thanks to a dermatographic pencil (when the transillumination technique is not available, the treatment area can be estimated at a 3mm distance posterior to the limbus).

Laser Settings

Recommended settings: Quantel SubCyclo™ mode duty cycle at 31.3% (0.5ms on, 1.1ms off) or duty cycle at 25% (0.62ms on, 1.9ms off), at 2000mW.

Placement (Fig. 4)

Unscrew the footplate carefully from the probe exposing only the metal threaded rod and ball fiber.

Put the tip of the probe in contact with the sclera and place it over the treatment area previously identified.

Hold the probe perpendicular to the surface of the eye.

Applications and Treatment (Fig.5)

To treat the superior hemisphere, apply firm pressure and move the device in a continuous sliding motion from the 9:30 to 2:30 clock positions (laser application area around the Iris).

Apply laser continuously, for 4 to 7 complete passes between clock positions, for a total of 50 to 80 seconds.

The 3 and 9 o'clock meridians must imperatively be avoided.

To treat the inferior hemisphere, repeat process for the suggested time over 3:30 to 8:30 clock hours.

Total treatment duration: 50 to 80s by hemisphere, equivalent to 100 to 160s by eye.

Performances, mode of action and clinical benefits

The silica fiber of the SUBCYCLO® devices allows the transmission of a laser beam from a vitrectomy machine to the operation site. Their specific transparent tip can be put against the patient's sclera to help stabilizing the laser beam. The laser beam is applied to the ciliary bodies, through the sclera, to decrease the aqueous humor production and, consequently, decrease the intraocular pressure.

No clinical performances or clinical benefits are associated with the use of SUBCYCLO®.

Compatibility with other devices

SUBCYCLO® is compatible with Quantel Medical Vitra 810 and Supra 810 lasers.

Information to be communicated to the patient

Any serious incident that occurs in relation with the device should be notified to FCI S.A.S. and the competent authority of the State in which the patient is established.

The user must inform the patient of the contraindications and potential complications related to the device that are set out in the present document.

Waste disposal

With household waste or recycle	With hospital waste (contaminated products)
Outer box, instructions for use	SUBCYCLO®
Pouches	

The disposal of this device does not present any physical risk.

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